

Improved T Cell Activation Bioassays for Development of Bispecific Antibodies and Engineered T Cell Immunotherapies

Pete Stecha, Denise Garvin, Jim Hartnett, Frank Fan, Mei Cong and Zhi-jie Jey Cheng

Promega Corporation, 2800 Woods Hollow Rd, Madison, WI 53711




1. Introduction

Immunotherapy aims to boost a patient's own immune system to fight disease. In recent years, a variety of immunotherapy strategies aimed at inducing, strengthening or engineering T cell responses have emerged as promising approaches for the treatment of cancer and autoimmune disease.


Here we describe a platform of T cell activation bioassays for the development of CD3 bispecific antibodies and engineered T cell immunotherapies.

Specifically, we developed two bioluminescent reporter-based bioassays to measure T cell activation via CD3 (NFAT-RE) or CD3 + CD28 (IL-2 promoter). These bioassays include the following:

TCR/CD3 (NFAT) effector cells: Jurkat cells engineered with an NFAT-RE driving luciferase expression. Responds to CD3, but not CD28 stimulation.




Lineberry-N et al. Nature Reviews Immunology 2007




TCR/CD3 (IL-2) effector cells: Jurkat cells engineered with an IL-2 promoter driving luciferase expression. Responds to CD3 and CD3+CD28 stimulation.

4. Abatacept Modulates CD3+CD28 T Cell Activation using TCR/CD3 (IL-2) Cells



Experimental Design


- (1) TCR/CD3 (NFAT or IL-2) effector cells are incubated with Raji (CD80/86+) target cells
- (2) T cell activation is induced via crosslinked anti-CD3 Ab and CD28 engagement by its ligand CD80/86 expressed on the Raji target cells
- (3) Addition of a CTLA-4/IgG fusion protein (Abatacept) binds CD80/86 and inhibits CD28-mediated T cell activation.



Increasing concentrations of Abatacept were added to either TCR/CD3 (NFAT) or TCR/CD3 (IL-2) effector cells, as indicated. Abatacept induced a significant decrease in TCR-mediated luciferase activity in TCR/CD3 (IL-2) effector cells compared to TCR/CD3 (NFAT) effector cells. This is expected because CD28 functions independently of the NFAT response element (see Introduction).

5. Measurement of CD3 Bispecific Antibody Activity: Assay Design and Protocol


Assay Design for Measuring CD3 Bispecific Antibody Activity



Experimental Design


- (1) TCR/CD3 (NFAT or IL-2) effector cells are incubated with increasing concentrations of a CD3 bispecific Ab
- (2) The bispecific Ab simultaneously binds to TCR/CD3 on the effector cells and tumor antigen on the target cells
- (3) Bispecific Ab binding stimulates IL-2 or NFAT luciferase activity

Assay Protocol for Measuring CD3 Bispecific Antibody Activity




2. Both TCR/CD3 (NFAT) and TCR/CD3 (IL-2) Cells Respond to TCR/CD3 Stimulation

T Cell Activation Protocol




TCR/CD3 (NFAT)

EC50=82ng/ml
Fold Induction=165




TCR/CD3 (IL-2)

EC50=35ng/ml
Fold Induction=76



TCR/CD3 (NFAT) (Left) and TCR/CD3 (IL-2) (Right) effector cells were stimulated with increasing concentrations of an anti-CD3 Ab.


3. T Cell CD3+CD28 Co-stimulation Measured using TCR/CD3 (IL-2) Effector Cells




■ anti-CD3 + anti-CD28 Abs: EC50=36ng/ml, Fold Induction=230
■ anti-CD3 Ab: EC50=11ng/ml, Fold Induction=31
■ anti-CD28 Ab: EC50=>368ng/ml, Fold Induction=8

TCR/CD3 (IL-2) effector cells were stimulated with increasing concentrations of anti-CD28, anti-CD3 or a combination of anti-CD3+anti-CD28 Abs, as indicated.

6. Analysis of Blinatumomab CD3 x CD19 Bispecific Antibody Activity



Blinatumomab belongs to a class of constructed monoclonal antibodies, bi-specific T-cell engagers (BiTEs), that exert action selectively and direct the human immune system to act against tumor cells. Blinatumomab specifically targets the CD19 antigen present on B cells.



Increasing concentrations of Blinatumomab were added to either TCR/CD3 (IL-2) or TCR/CD3 (NFAT) effector cells, as indicated. Blinatumomab induced a dose-dependent increase in luciferase in both TCR/CD3 (IL-2) and TCR/CD3 (NFAT) effector cells in the presence of Raji (CD19+) target cells. No response was detected in the absence of Raji (CD19+) target cells.

7. Assay Qualification with Blinatumomab: Assay Precision, Accuracy and Linearity

Accuracy and Intermediate Precision (N = 6)


Assay Qualification Design:

- Two analysts
Three days
Four plates per day
- 100% vs 50%
 - 100% vs 75%
 - 100% vs 150%
 - 100% vs 200%


Data shown are generated using TCR/CD3 (NFAT) effector cells, blinatumomab, and Raji (CD19+) target cells.

Expected Relative Potency	Assay	Analyst	Measured Relative Potency %	Mean %	SD%	Accuracy Recovery %	Precision n	RSD %
50%	1	1	45.3	47.8	3.6	95.7	7.6	
	2	1	45.0					
	3	1	45.5					
	4	2	46.2					
	5	2	52.2					
	6	2	52.8					
75%	1	1	62.2	71.2	6.8	95.0	9.5	
	2	1	73.9					
	3	1	63.3					
	4	2	78.4					
	5	2	75.9					
	6	2	73.6					
150%	1	1	144.5	134.4	15.5	89.6	11.5	
	2	1	143.5					
	3	1	121.4					
	4	2	108.8					
	5	2	143.4					
	6	2	144.7					
200%	1	1	174.9	184.8	15.0	92.4	8.1	
	2	1	195.5					
	3	1	179.8					
	4	2	162.0					
	5	2	198.9					
	6	2	198.0					
Overall							93.2	9.2

Repeatability (%CV) = 3.01%




Linearity and Range (R²=0.995)




8. Measurement of Chimeric Antigen Receptor T (CAR-T) Cell Activity

CAR-T Design




TCR/CD3 (NFAT)



Experimental Design

- (1) TCR/CD3 (NFAT or IL-2) effector cells were engineered to express anti-CD19/CD20 CAR-T receptors
- (2) The effector cells were incubated in presence or absence of Raji (CD19/CD20+) target cells

TCR/CD3 (IL-2)



Serial dilutions of anti-CD19/CD20 CAR-T effector cells were incubated in the presence and absence of Raji (CD19/CD20+) target cells. Luciferase activity was detected in the presence of Raji cells, but not with the CAR-T effector cells alone.

9. Conclusions

We have developed a platform of T cell activation bioassays that incorporate a bioluminescent reporter-based readout of T cell activation via CD3 (NFAT-RE) or CD3 + CD28 (IL-2 promoter). These assays reflect the mechanisms of action of biologics designed to engage, recruit, and stimulate T cell activation to attack target disease cells. Specific applications include measurement of **anti-CD3 bispecific Ab** and **CAR-T cell activity**.

The bioassays provide the following:

- Mechanism of action (MOA)-based measure of biologics activity**
 - Specific measure of CD3 or CD3 + CD28 T cell activation pathways
 - Quantitative measure of anti-CD3 Ab and bispecific Ab potency
- Consistent and reliable measure of biologics activity**
 - Demonstrated precision, accuracy, reproducibility, robustness
 - All assays can be used as "Thaw-and-use" cell format, no cell culture required
 - Functional performance suitable for development into potency, stability, and NAB assays (data not shown)
- Easy-to-implement**
 - Rapid and convenient workflow
 - Amenable to standard 96-well and 384-well plate formats