

Laboratory Refrigerators B Medical Systems | L Range

Laboratory Refrigerators are devices intended for the safe storage of samples, specimens, cultures and other laboratory preparations at cold temperatures between +4°C and +15°C.

Compliant to DIN 13221 | Medical Device Class I according to MDD 93/42/EEC



medical systems

SAVING LIVES THROUGH RELIABLE AND INNOVATIVE TECHNOLOGY

Medical Refrigeration



Safety Standards B Medical Systems



The Safety Standards developed by B Medical Systems define certain significant technical features of a product. These ensure the safe storage of the preparations as well as setting the highest standards of safety for the user.

PRECISION LINE	PREMIUM LINE	
		B Medical Systems Electronics
		B Medical Systems digital display
		B Medical Systems 4,3" digital color display with capacitive buttons
	0	B Medical Systems 7" full touchscreen display
		Safety door lock (with 2 keys) and key-operated power switch ON/OFF (with 2 keys)
		Power indicator light and digital temperature indicator (display: 0.1 digits)
		Controlled fan cooling system for constant temperature and even temperature distribution across the entire refrigerating chamber [not applicable for model LF260 (Freezer compartment)]
		Automatic switch-off of the evaporator fans when the door opens [not applicable for model LF260 (Freezer compartment)]
		Self-contained alarm system with integrated battery takes over the alarm function and temperature value measurements in case of system failure for at least 48 hours
		Acoustic and visual alarm signal in case of temperature alarm and system failure
		The alarm history function on the electronic stores all the relevant values during a temperature alarm, such as: min., max., average temperature and alarm duration
		Remote transmission alarm signal (via potential-free contact) in case of temperature alarm (change-over contact)
		Safety thermostat prevents dropping of the temperature below +2°C
		Door opening alarm (visual / acoustic)
		Designed and tested for climatic class SN (ambient temperature range +10°C to +32°C)
		Designed and tested for climatic class SN / T (ambient temperature range +10°C to +43°C)
		Interior made from stainless steel
		Internal LED lighting
		Additional remote transmission alarm signal (via potential-free contact) in case of system failure (change-over contact)
		Smooth castors for optimum flexibility of movement
	0	Ambient temperature sensor
\bigcirc		°B Connected - Universal software for the monitoring of refrigeration devices, including the acquisition, recording and visualization of temperature data
\bigcirc		DCU (Data Communication Unit) - Hardware module monitoring all operating conditions and passing them through to a central database – via local network (in combination with °B Connected)
		Ethernet interface for the visualization of all operating and control functions (hardware and software settings) via °B Connected monitoring software on a peripheral device (computer)

MEDICAL REFRIGERATION L RANGE | Laboratory Refrige



Laboratory Refrigerators B Medical Systems | L Range

Laboratory Refrigerators are devices intended for the safe storage of samples, specimens, cultures and other laboratory preparations at cold temperatures **between +4°C and +15°C** with multiple safety features developed as a result of over 35 years of expertise in refrigeration technology.

9 models • Volume 52 > 895 L • Set temperature +5°C • Climate class SN | SN/T • Compliant to DIN 13221 | MDD 93/42/EEC, Class I

In conformity with national and international guidelines, regulations for Medical Devices offering reliability, efficiency and safety at an optimal price.



• Integrated multi-function monitoring electronics

Integrated B Medical Systems multi-function electronics; with easy access in the door level; offers alarms, central alarm system, set point security, battery backup, full display and many other control features; compatible with °B Connected monitored solution for data monitoring and recording.

• Lifetime warranty corrosion-free cabinet

• One piece material cabinet with a lifetime warranty for zero corrosion, highest quality polypropylene for better insulation and no sharp edges to provide easier maintenance and cleaning

 Smart ergonomy - Easier and safer handling because of heavier components placed at the bottom

• High storage capacity

High storage per square meter in terms of net volume, with a very versatile inner volume utility due to multiple user friendly rails allowing modularity.



MODELS L55-130-290-380 & LF260 | PRECISION LINE



• Exclusive integrated electronics handle bar (Patent pending)

All functionalities are easily accessible:

- Digital control with full functionalities -Standard multifunction digital sensitive 4.3" panel integrated at optimal level in the door handle with full commands and connections, or optional full 7" touchscreen with pre-installed connection allowing exclusive °B Connected monitoring functionalities
- Open/close Heavy duty door lock mechanism

• High quality materials

New high quality coating, certified medical devices quality and antibacterial, high quality steel for better longevity and easy hygiene control.

Versatile modularity

Excellent storage capacity and modularity - Large choice of shelves and drawers offering modularity for every need.





MEDICAL REFRIGERATION L RANGE | Laboratory Refrige





	L	55	Ľ	130	Lä	290	L3	880	LF	260	Refrigerator 8 combined	Freezer
Gross / Net volume (I)	52 / 45		121/106		297 / 269		362 / 329		R 121 / 1 0)6	F 121 / 10	6
Set temperature (preset)					+	-5°C					-32°C	
Set temperature (setting range) can be adjusted in steps of 0.1°C					+4°C	to +15°C					-32°C to -2	0°C
Temperature cold / warm alarm limit					+3°(C / +7°C					-37°C / -27	°C
Hold over time (+5°C to +10°C)	1 h 15		1 h 27		59 min.		54 min.		1 h 27		42 min.	
Climate class (ambient temperature range)					SN (+10°	C to +32°C)						
Defrosting technique					Na	itural					Manual (for	freezer)
Refrigerant type					R	500a					R290	
External dimensions H x W x D (mm)	670 x 495 x	< 575	830 x 595 x	x 695	1670 x 595	x 695	1980 x 595	x 695	1660 x 595	x 695		
Inner dimensions H x W x D (mm)	470 x 375 x	< 350	630 x 475 x	k 470	1470 x 475	x 470	1790 x 475	x 470	630 x 475 >	470	630 x 475 >	470
Net weight with standard equipment (kg)	39		59		94		107		117			
Supply voltage (V)	220-240	115-127	220-240	115-127	220-240	115-127	220-240	115-127	220-240	115-127	220-240	115-127
Frequency (Hz)	50 / 60	60	50 / 60	60	50 / 60	60	50 / 60	60	50 / 60	60	50 / 60	60
Power (W)	100	100	200	200	200	200	200	200	200	200	300	360
Gage Energy consumption (kWh/24h)	0.38 / 0.45	0.38	0.89 / 1.13	0.61	0.96 / 1.27	1.06	0.93 / 1.26	1.10	0.89 / 1.13	0.61	1.35 / 1.28	1.98
Heat emission (Kcal/h)	9	8	29	18	29	22	37	20	29	18	53	65
Compressor running time (%)	25	20	31	24	28	30	36	25	31	24	37	36
Noise level (dB(A)) (at 1m height & 1m distance)	35 / 36	37	35 / 36	38	38 / 40	39	40	38	35 / 36	38	45 / 47	47



		L	400	LS	500	L7	00	LS	00		
Gro	oss / Net volume (I)	451/384		598 / 510		747 / 636		895 / 763			
Set	temperature (preset)					+5°C					
	temperature (setting range) be adjusted in steps of 0.1°C		+4°C to +15°C								
Ten	nperature cold / warm alarm limit	+3°C / +7°C									
Но	ld over time (+5°C to +10°C)	1 h 36		1 h 36		1 h 36		1 h 36			
Clir	nate class (ambient temperature range)	SN / T (+10°C to +43°C)									
Dej	frosting technique				N	atural					
Ref	rigerant type				F	R600a					
Fxt	ernal dimensions H x W x D (mm)	1988 x 699	x 1039	1988 x 845	x 1039	1988 x 992	x 1039	1988 x 113	9 x 1039		
	er dimensions H x W x D (mm)	1167 x 447		1167 x 590	x 713	1167 x 740	x 713	1167 x 887	x 713		
	t weight with standard equipment (kg)	194		215		237		259			
	Supply voltage (V)	220-240	115-127	220-240	115-127	220-240	115-127	220-240	115-127		
	Frequency (Hz)	50 / 60	60	50 / 60	60	50 / 60	60	50 / 60	60		
y Save	Power (W)	300	300	300	300	300	300	300	300		
Energy	Energy consumption (kWh/24h)	1.40 / 1.40	1.30	1.40 / 1.40	1.30	1.40 / 1.50	1.30	1.50 / 1.50	1.30		
	Heat emission (Kcal/h)	29	27	33	28	32	30	35	29		
	Compressor running time (%)	22	22	25	23	24	25	25	24		
	Noise level (dB(A)) (at 1m height & 1m distance)	46 / 49	52	47 / 51	50	50/51	46	47 / 47	46		

rators

All values were measured at +25°C ambient temperature and without load (with inertial mass). Subject to change without prior notice. Some of the accessories shown in the pictures are optional.

.)

REFRIGERATION	Laboratory Refrigerators
MEDICAL I	L RANGE

RANGE

Technical Data Specifications





B

		L55-130 & LF260	L290-380	L400-500-700-900						
Regulation se Precision (from -80	n sor 0°C to +180°C) ± 0.3°C	NTC		PT1000						
Display sensor Precision (from -80°C to +180°C) ± 0.3°C in reference body / 100 ml DOW corning 200-5CST (Silicon Oil)		NTC		PT1000						
Accu data / function time of the control panel when system failure		12 V - 2.3 AH	/ 48h	12 V - 7 AH / 48h						
Relative humidity at +32°C		≤ 75%								
Door insulatio	on (polyurethane)	80 mm P	J	65-82 mm PU						
Casing insulation (polyurethane)		60 mm P	J	77-82 mm PU						
	Inner body	Polypropylene (cop	olymer UV)	Stainless steel (1.4301)						
	Outer body	Polypropylene (cop	olymer UV)	DX51D + Z100 coated						
Matarial	Inner door	Polypropylene (copolymer UV)	Aluminium	Stainless steel (1.4301)						
Material	Outer door	Polypropylene (cop	olymer UV)	DX51D + Z100 coated						
	Drawer	Polycarbonate, tra	ansparent	Stainless steel (1.4301)						
	Shelf		Wire DIN1	77, PA11 coated						
European Me	dical Device Directive		93 / 42 / EEC, Cla	SS (models: 220-240 V - 50/60 Hz)						
EMC directive	?		2014 / 30 / EU	(models: 220-240 V - 50/60 Hz)						
Low voltage a	lirective	2014 / 35 / EU (models: 220-240 V - 50/60 Hz)								

Yes (models: 115-127 V - 60 Hz)

Class ISO 6 / EC GMP B

FDA regulation Class I medical devices acc. to 21CFR Part 862.2050 (models: 115-127 V - 60 Hz)

cCSAus certificate

GMP - Clean room classification



Equipment Standard & optional

RANGE

		LSS	L130	L290	L380	LF260	L400	L500	L700	L900
B Medical	Digital display	•	•	•	•	• 2	-	-	-	-
Systems	4.3" digital color display	-	-	-	-	-	•	•	•	•
Electronics	7" full touchscreen display	-	-	-	-	-	0	0	0	0
	Drawer with Front Cover	01	01	()≤5	()≤6	02 01	-	-	-	-
	Drawer without Front Cover	() ≤ 3	() ≤ 3	()≤12	()≤15	() ≤ 3	-	-	-	-
Interior equipment	Bottom Drawer	•1	• 1	• 1	• 1	• 2	-	-	-	-
lyapment	ST-Drawer with Front Cover	-	-	-	-	-	()≤3	()≤3	()≤3	() ≤ 3
	Wire Shelf	• 1	• 2	• 5	6	• 2	• 5	• 5	• 5	• 5
Internal LED li	ghting	•	•	•	•	•	•	•	•	•
°B Connected	- Monitoring Software	0	0	0	0	0	•	•	•	•
DCU - Data Co	pmmunication Unit th °B Connected)	0	0	0	0	0	-	-	-	-
Ethernet inter	face	-	-	-	-	-	•	•	•	•
Ambient temp	perature sensor	-	-	-	-	-	0	0	0	0
Potential-free n case of system fo	alarm contact ^{ailure}	•	•	•	•	•	•	•	•	•
Integrated pol	rt for external sensor ^{mer)}	•	•	•	•	• 2	•	•	•	•
Rollers Smoo	oth castors with stabilizers	- -	• 2 -	• 2 -	• 2 -	• 2 -	- • 4	- • 4	- • 4	- • 4
Door hinge Rig	ght Left				$\bullet \mid \bigcirc$				$\bullet \mid \bigcirc$	



Accessories Standard & optional

	L55	L130	L290	L380	LF260	L400	L500	L700	L900
Dividers Set (H 35 mm or H 100 mm) to facilitate the handling and storage management	01	01	()≤5	()≤6	() ≤ 2	-	-	-	-
S-Rack (H 100 × W 113 × D 638 mm)	-	-	-	-	-	0	0	0	0
Temperature chart recorder pre mounted (in the form of a circular chart recorder)	0	0	0	0	0	0	0	0	0



Full functionalities at a glance

B Medical Systems | Electronics

PRECISION LINE | Multifunctional electronics with digital display and easy access, in the door handle; compatible with °B Connected monitored solution.




• Key power switch O = OFF / I = ON

- To switch the unit on, turn the key switch to position "I". The green power LED lights up



- After the self-check, the temperature inside the device is displayed
 - O Digital display
- During normal operation the display shows:
- Left-hand side: current time and internal temperature displayed in °C (Celsius) or °F (Fahrenheit) in increments of 0.1°
 - Right-hand side: the maximum and minimum inside temperature



• Acoustic and visual alarm signal

When an alarm occurs, the red alarm LED lights up and an acoustic signal sounds. The corresponding alarm message appears on the display alternately with the measured inside temperature. If several alarms occur simultaneously, the messages will be shown alternately.



O Function buttons

- UP / DOWN buttons: navigation within the parameters and modification of their values
- ENTER button: use the ENTER-key to confirm the actual selection. This applies to the parameter value. The change of a parameter only becomes effective when confirmed
- MUTE / BACK button: in case of an alarm, the acoustic alarm can be muted for a specified time. In Parameter-mode the same key acts a "BACK-button" to go one step back
- LIGHT button: allows to switch on or off the interior light of the device without opening the door

O USB-port

A USB port is available for data export only. Measured values of the display sensor are stored for at least 30 days (ring buffer) together with a time stamp.



/**+**_/

~

← 1(×

*

The new B Medical Systems Electronics

offers a wide range of adjustment and diagnostic facilities as well as additional protection / warning operations (via external alarm operations, histories and individual display signals).

PREMIUM LINE | Multi-function electronics with 4,3" digital color display and capacitive buttons, integrated at optimal level in the door handle with full commands and connections.



• Power switch OFF / ON

• To switch on press the ON/OFF button for 1 second. The green LED will light up

 Switch off protected by password. The door will automatically unlock before powering down

B cond 200 Tree 5.0°C solution 5.0°C solution

• Digital color 4,3" display

During normal operation the display shows:

• At the top: current voltage and current time

 Left-hand side: internal temperature displayed in °C (Celsius) or °F (Fahrenheit) in increments of 0.1°

• Right-hand side: temperature of the regulation sensor. Displaying this temperature can be disabled through the settings



OUSB-port & SD-card-slot



Discover the ultimate control experience with the **new B Medical Systems 7" full touchscreen display**

More comfortable, more functionalities, no compromise!

• Acoustic and visual alarm signal

When an alarm occurs, the red alarm LED lights up and an acoustic signal sounds. The corresponding alarm message appears on the display together with the measured inside temperature. If several alarms occur simultaneously, the messages will be shown alternately.

• Capacitive buttons

- ARROW buttons: select screen items with the arrow buttons. In parameter editing mode, the arrow buttons are used to change values
- ENTER button: activates the currently selected item
- BACK button: goes one level upwards in the menu tree structure. It also serves as "escape" button
- DOOR UNLOCK button: unlocks the door of the device
- MENU button: opens the menu screen
- ALARM MUTE button: deactivates the acoustic alarm for a predefined period of time



←

G

:=

II ×

Optimum control and protection

B Medical Systems | °B Connected & DCU

°B CONNECTED | MONITORING SOFTWARE

Universal software for the monitoring of refrigeration devices, including the acquisition, recording and visualization of temperature data.

- Unique monitoring software for the full range of Blood Management Solutions and Medical Refrigeration products
- Web-based interfaces for computers and mobile devices
- Modern design for simple and intuitive use
- Graphic display of temperature curves
- Integrated event and activity history of appliances' components
- Data recording on centralized database for long-term archiving
- Easy setting of specific alarm forwarding, via email or SMS
- Generation of reports compiling crucial data and activities
- Temperature and detailed device data export for third-party software
- Important price advantage compared to a traditional circular chart recorder and its spare parts
- REST API to access raw data directly from the database in read only mode

KEY BENEFITS:

- Unique interface for the monitoring of the full range of refrigeration devices
- Centralized database providing data access to entire customer network

	URATION	
B con	nected	+O USER MANAGEMENT
	ER MMI (HIGDNO)	 Password login Access right
ne : 1 Device configura	ation	+
PREMIUM LINE		DEVICE SETTINGS
		 Device information Device & user configuration
A AMERAN A AMERAN A AMERAN A AMERAN A AMERAN A AMERAN		((ر٢٠))
CST ed. (m.e.)	88 mm 2	ALARM PROFILE
BATOR		 Parameters of the alarm profile (day, hour, minute) Alarm management

2 DISPLAYING




DETAILED DISPLAY

 Optimised design for mobile and touch compatible devices
 Support of multiple languages

 Real time data: Current temperatures, set point, upper and lower limit, door state, compressor runtime
 Log data: Events, alarms within date/time, event type, description of the event. Event and alarm confirmation function
 Batch management: Read in of plasma bags with a barcode scanner
 Movement state and door control of Agitator devices
 Parameter settings of each device: Factory settings, real time clock, language, firmware version. R5485 address. IP address °B Connected & DCU offer a highly flexible system that is adaptable to specific customer requirements

- > Complete & legally safe documentation of temperature data
- > Comprehensive applications and diagnostic possibilities



DCU DATA COMMUNICATION UNIT (PRECISION LINE)

Hardware module monitoring all operating conditions and passing them through to a central database – via local network.



- Interface connection of B Medical Systems appliances to an existing network
- Graphical user interface displaying temperature and configurating devices and alarms easily
- Direct connection to Ethernet and serial bus RS485
- Digital in/out (programmable and customer-specific use)
- Recording and storing appliance relevant data
- Integrated USB port allowing to export stored data
- The DCU combined with the °B Connected software replaces the paper temperature recorder
- All data are recorded and save in the internal storage of the DCU and backed up in °B Connected data base if connected
- Several additional self-sufficient temperature sensors (up to 4 PT1000) may be connected
- Humidity sensor input (4-20 mA)

KEY BENEFITS:

- Central system collecting relevant temperature data of the appliances and their respective operating conditions
- A number of connection facilities allowing flexible upgrades for individual projects



*

Safe global blood management: from collection to transfusion, transportation, processing and storage



Reliable solutions for safe vaccination around the world



State-of-the-art technology for the exacting needs of the medical world

Our Global Expertise



After Sales support and service

We strive to provide you with the highest standards of service; not only through our selected distributors and partners for all your maintenance and service but also our second line trouble shooting and after sales service. This Factory based group of engineers is there to help our partners and yourself to get the best solution for your cold storage needs.



SAVING LIVES THROUGH RELIABLE AND INNOVATIVE TECHNOLOGY

B Medical Systems (formerly Dometic Medical Systems) has more than 35 years' experience in the medical refrigeration sector.

The company, formerly known as Electrolux Medical Systems, was founded in 1979 when the World Health Organization approached Electrolux in Vianden, Luxembourg, to create a solution for the safe storage and transport of vaccines around the world. In 2001, Electrolux Medical Systems became part of the Dometic Group, and was renamed Dometic Medical Systems. Having established a legitimate reputation in the medical equipment industry, the company has also become a global leader in vaccine cold chain.

B Medical Systems S.à r.l.

17, op der Hei L - 9809 Hosingen, Luxembourg

Tel.: (+352) 92 07 31-1 Fax: (+352) 92 07 31-300 info@bmedicalsystems.com









Ed 09/2018

B MEDICAL SYSTEMS; the °B logo; Ultra SDD and Saving lives through reliable and innovative technology; are trademarks of B Medical Systems S.à r.l. | B Medical Systems products are protected by granted and pending patents and registered designs. All rights reserved. All indications without guarantee and excluding all liability due to faulty, incomplete or outdated information. All technical data subject to change without prior notice.