

RADSIGHT Access

PART NUMBER: 43700100013


INTENDED USE

SMALL HOSPITALS, SMALL NUCLEAR FACILITIES, OR AS A TRIAL OR STARTER SYSTEM BEFORE UPGRADING TO A LARGER SYSTEM.

Package to provide a dosimetry solution as a software only package with an entry level price. This is an ideal system for getting started or for small scale use. Can be installed on a single computer for stand alone operation or as a server installation with a separate computer running the RadSight application to issue and return EPDs. Dosimetry data can be viewed from any computer on the network using a web browser.

WHAT'S INCLUDED

- RadSight Access software licensed for up to 10 EPDs
- Installers for RadSight server component
- Ability to add 1 issue/return station



EXAMPLE

Single computer installation for simple dosimetry.

Installation requirements:

RadSight Windows 10 can also be installed on any Windows 10 computer as long as it meets the minimum versions for Windows 10 shown below.

RadSight Windows 7 installation requirements:

- Windows 7 Service Pack1 (x64)
- Windows Server 2016 (minimum version 1709)
- Windows 10 Anniversary Update (version 1607) (x64)
- Windows 10 Creators Update (version 1703) (x64)
- Windows 10 Fall Creators Update (version 1709) (x64)
- Note that 32 bit versions of Windows 7 or Windows XP are not supported

Thermo Fisher Scientific is an authorized reseller of Remote DNA, RadSight Access and RadSight Live.

For more information, contact your Thermo Fisher sales representative or visit us online.

Phone: 1-800-274-4212 (NA); +49 9131 998 171 (outside NA) | Email: customerservice.rmsi@thermofisher.com | www.thermofisher.com/epdtrudose

BASE