



**Thermo Scientific**  
Chromatography Columns  
and Consumables 2016-2017


# Connected chromatography solutions

**Thermo**  
SCIENTIFIC

# Sample Preparation Products

Your sample is precious and you cannot afford failures. We have a comprehensive range of products delivering clean and consistent extractions using innovative sample preparation technologies for your analytical needs.

## Featured Products



### SMART Digest Kits

Designed for biopharmaceutical and proteomic applications that require highly reproducible, sensitive and fast analyses of proteins often in high throughput workflows.

**PAGE 1-004**



### SOLA Products

This first fritless SPE product range which has been specifically designed for bioanalysis provides greater reproducibility with cleaner, more consistent extracts. Thermo Scientific™ SOLA™ products provide unparalleled performance characteristics compared to conventional SPE, phospholipid removal and protein precipitation products.

**PAGE 1-007**



## Section Contents

SMART Digest Kits . . . . .	<b>1-004</b>	HyperSep SPE Cartridges and Plates. . . . .	<b>1-016</b>
SOLA SPE Products . . . . .	<b>1-007</b>	SPE Manifolds and Accessories . . . . .	<b>1-025</b>
HyperSep SLE Cartridges and Plates. . . . .	<b>1-012</b>	HyperSep Dispersive SPE (QuEChERS) . . . . .	<b>1-030</b>
Microscale SPE Plates and Tips. . . . .	<b>1-014</b>	Syringes and Filters . . . . .	<b>1-034</b>



### QuEChERS

A quick and easy sample preparation approach for pesticide analysis. Convenient extraction and clean-up of pesticides in food and other complex matrices.

**PAGE 1-030**



### Syringe Filters

Providing cleaner sample extracts by removing interfering materials and fine particles. Delivering high quality, robust solutions for a wide range of samples and applications.

**PAGE 1-034**

# Sample Preparation

## Technical Resources



visit [www.thermoscientific.com/chromexpert](http://www.thermoscientific.com/chromexpert)  
to access the following information:

- SPE Phase Selection
- SPE Phase Selection by Manufacturer
- SPE Procedure - Six Steps for a Clean Extract
- Method Development Optimization in SPE



# SMART Digest Kits

## Facilitating perfect digestion

Thermo Scientific™ SMART Digest™ kits are designed for biopharmaceutical and proteomic applications that require:

- High analytical reproducibility
- High sensitivity
- Fast analysis
- High throughput workflows

They achieve this due to the optimized heat stable immobilized trypsin design. It offers significant benefits over existing conventional in-solution digestion protocols:

- Simple process
- Significantly faster than in-solution digestion
- More reproducible than in-solution digestion
- Higher sensitivity
- Higher throughput
- Ability to automate




### Simplicity of use with SMART Digest Kits


Current sample preparation protocols for the digestion of proteins are multifaceted and laborious. Without a single source solution, errors are far more likely to occur. Therefore, the overall protocol is subject to irreproducibility and significant sample processing time.

The SMART Digest kit overcomes these barriers by providing a simple to implement, integrated process (see the diagram below), which provides high data confidence and significantly increased reproducibility. The entire protocol can be easily automated for high throughput processing.

### Standard Digest




### SMART Digest



### Fast digestion with SMART Digest Kits

The SMART Digest process is not only simple, but significantly reduces the time from sample preparation to analysis. Typically it takes less than 60 minutes, dependent on sample complexity; to achieve full digestion as can be seen in the IgG example below, full digestion is achieved in 60 minutes.

#### IgG Assay m/z- 937




#### Recommended digestion starting conditions for known proteins\*

Protein	Digest Time (min)
Insulin	4
BSA	< 5
Carbonic anhydrase	< 5
Lysozyme	< 5
Apo-B	30
IgG	45
IgG in 50 µL plasma	75
Ribonuclease A	150
Thyroglobulin	240
C-reactive protein	240

\*200µL protein solutions were used for all samples. A concentration of 100µg/mL was used for all samples, except IgG in plasma which used 17.5mg/mL

### SMART Digest vs. conventional in-solution digest protocol

#### In-solution digest



#### SMART Digest




- Sample preparation
- Protein assay
- Digestion
- SPE + dry/reconstitute
- Spin down

### Reproducibility of digestion with the SMART Digest Kits

The SMART Digest kit provides significant improvements in reproducibility over existing protocols, which results in fewer sample failures, higher throughput and the ability to more easily interrogate data.

The chromatogram to the right shows two overlaid traces of BSA digests produced by two separate users who had never performed a digestion protocol before.




### Higher sensitivity with the SMART Digest Kits

The SMART Digest kit can provide higher sequence coverage due to more efficient digestion and the reduced number of chemicals used, resulting in less ion suppression during LC-MS analysis. The reduced use of chemicals also decreases the number of chemically-induced post-translational modifications.

The mass spectrometry (MS) peak area response comparison on the right shows how Thyroglobulin, digested using the SMART Digest kit, provides higher levels of sensitivity compared to a standard in-solution digest. In addition to this the SMART Digest protocol took two hours compared to twenty hours for the in-solution protocol.

### Comparison to Historical Method



\*Clarke, N. J., et al.; *J. Inv. Med.* **2012**, *60*, 1157-1163

### SMART Digest Kits

Description	Cat No.
Smart Digest kit	<b>60109-101</b>
Smart Digest kit with 96 well filter plate	<b>60109-102</b>
Smart Digest kit with SOLA $\mu$ HRP SPE plate	<b>60109-103</b>

Find out more about how SMART Digest Kits can benefit your workflow at:  
[www.thermoscientific.com/SMARTdigest](http://www.thermoscientific.com/SMARTdigest)



# SOLA SPE Products

## Consistent excellence for bioanalysis


Thermo Scientific SOLA products are designed for bioanalytical and clinical research analysts who are tasked with providing high quality analytical results from complex biological samples in a high throughput environment, while complying with strict legislation. These demands are compounded by the continued push to higher efficacy drugs and long acting formulations which continue to drive sensitivity requirements to lower levels to enable accurate quantification.



In order to meet these demands, bioanalytical methods must provide:

- Robustness – low analytical failure rates
- Ability to process low sample volumes
- High sensitivity
- High reproducibility
- Ease of use
- High throughput processing
- Efficient and fast processes

The SOLA SPE range meets these demands due to the unique and innovative frit-less SPE technology which eliminates the issues with traditional loose-packed SPE formats (see figure to the right). Combining the support material and active media components into a solid, uniform sorbent bed provides stable and controllable flow through characteristics and has an added advantage when dealing with viscous biological samples, as it prevents blocking and enables high throughput processing.




SOLAμ SPE design – limiting issues associated with conventional SPE formats

The manufacturing process also allows for high levels of reproducibility, not only from cartridge to cartridge or well to well but also batch to batch (see figure to the right).

The SOLAμ SPE range has the added benefit of being able to provide:

- Up to a 20 fold increase in sensitivity
- Ability to process samples restricted in volume
- Increased workflow efficiency and sample integrity




Batch to batch reproducibility of SOLA compared to a conventional SPE product

### Increase your MS sensitivity


Up to 20 times increase in sensitivity can be achieved by loading a large volume of sample and eluting in a low volume.

In the following example 500µL of human plasma was loaded onto the SOLAµ plate for the analysis of niflumic acid. The compound was eluted in 25µL, providing a 20 times increase in concentration whilst maintaining excellent precision.



**The problem:**  
 improvement in assay sensitivity required

**SOLAµ solution:**  
 SOLAµ for up to 20 times increase in concentration of sample without changes to workflow







	Precision Data for Niflumic Acid Peak Area Ratio (%RSD) n = 18	Recovery of Niflumic Acid (%)	Matrix Effects (%)
QC Low (0.4ng/mL)	1.31	89.9	8.63
QC High (30ng/mL)	1.06	94.0	3.21

## Sample limited assays

SOLA $\mu$  allows you to directly scale down the volumes used in analytical methods, allowing for a reduction in sample usage and eliminating issues caused by evaporation, without compromising the sensitivity of your assay. This is also an important consideration when sample volumes are limited.

The following example shows that by loading 25 $\mu$ L of niflumic acid sample onto the SOLA $\mu$  plate and eluting in a total volume of 25 $\mu$ L, a ten-fold decrease in sample volume was achieved when compared to a traditional scale higher bed weight product. Equivalent method performance and high levels of reproducibility provided by SOLA technology were still maintained.



Equivalent peak response achieved for Niflumic Acid with SOLA $\mu$  using 25 $\mu$ L of sample compared with 250 $\mu$ L of sample with a 10mg SOLA SPE device


Precision Data for Niflumic Acid		
	Analyte Peak Area (%RSD)	Peak Area Ratio (%RSD)
Low QC	7.32	0.356
High QC	5.33	0.195

### Increased workflow efficiency


With traditional SPE the eluted sample is typically dried down to increase the concentration of the sample and thus improve the sensitivity. This causes an issue for certain compound types which can be lost during this step, resulting in reduced sensitivity.

SOLA $\mu$  allows the sample to be extracted without the need for dry down and reconstitution. Not only does this maximize recovery of the analytes it also improves workflow efficiency and increases productivity.


In the case of extraction of ibuprofen a four-fold pre-concentration was achieved without the need for dry down by loading 200 $\mu$ L of sample onto the SOLA $\mu$  plate and eluting in a total of 50 $\mu$ L. The results demonstrate that even with this low elution volume, excellent reproducibility was achieved.



**The problem:**  
workflow requires dry down to achieve sensitivity



**SOLA $\mu$  solution:**  
low elution volume allows removal of dry down



	Ibuprofen (%RSD) n=18	Ibuprofen recovery (%)
Low QC (25ng/mL)	4.00	90
High QC (750ng/mL)	1.70	95

### SOLA and SOLA $\mu$ meet bioanalytical needs by providing:

- A robust low sample volume preparation platform
- Reproducibility at low sample and solvent levels
- Increased workflow efficiency and throughput
- Sample enrichment (20 times)
- Mitigation against solvation and non-specific binding issues

**SOLA $\mu$  96 Well Plates**

Description	Bed Weight (mg)	Column Volume (mL)	Cat No.	Quantity
SOLA $\mu$ HRP	2	1	<b>60209-001</b>	1 Each
SOLA $\mu$ SCX	2	1	<b>60209-002</b>	1 Each
SOLA $\mu$ SAX	2	1	<b>60209-003</b>	1 Each
SOLA $\mu$ WCX	2	1	<b>60209-004</b>	1 Each
SOLA $\mu$ WAX	2	1	<b>60209-005</b>	1 Each

**SOLA SPE Cartridges**

Description	Bed Weight (mg)	Column Volume (mL)	Cat No.	Quantity
SOLA HRP	10	1	<b>60109-001</b>	100 Pack
SOLA SCX	10	1	<b>60109-002</b>	100 Pack
SOLA SAX	10	1	<b>60109-003</b>	100 Pack
SOLA WCX	10	1	<b>60109-004</b>	100 Pack
SOLA WAX	10	1	<b>60109-005</b>	100 Pack

**SOLA 96 Well Plates**

Description	Bed Weight (mg)	Well Volume (mL)	Cat No.	Quantity
SOLA HRP	10	2	<b>60309-001</b>	1 Each
SOLA SCX	10	2	<b>60309-002</b>	1 Each
SOLA SAX	10	2	<b>60309-003</b>	1 Each
SOLA WCX	10	2	<b>60309-004</b>	1 Each
SOLA WAX	10	2	<b>60309-005</b>	1 Each

# HyperSep SLE Cartridges and Plates

Solid supported liquid/liquid extraction (SLE) is a fast and effective sample preparation technique that provides considerable benefits over liquid-liquid extraction (LLE) protocols for removal of phospholipids from biological samples.

SLE offers the following advantages:

- Greater reproducibility and recoveries compared to LLE techniques
- Prevents emulsification often associated with LLE
- Reduced solvent requirements compared to LLE
- Can be completely automated unlike LLE
- Improved cleanliness of sample extract compared to protein precipitation techniques
- Improved sensitivity compared to protein precipitation techniques



Thermo Scientific™ HyperSep™ SLE is available in cartridge and 96 well plate formats and in acidified or basified formats to help simplify the protocol when acidic or basic mobile phases are used.

## HyperSep SLE Cartridges

Special treated diatomite SLE	Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
pH=7	200	3	<b>60109-200-3-7</b>	50 Pack
	500	3	<b>60109-500-3-7</b>	50 Pack
	500	6	<b>60109-500-6-7</b>	30 Pack
	1,000	6	<b>60109-1000-6-7</b>	30 Pack
	2,000	12	<b>60109-2000-12-7</b>	20 Pack
	4,000	25	<b>60109-4000-25-7</b>	15 Pack
	20,000	60	<b>60109-20000-60-7</b>	10 Pack
pH=9	200	3	<b>60109-200-3-9</b>	50 Pack
	500	3	<b>60109-500-3-9</b>	50 Pack
	500	6	<b>60109-500-6-9</b>	30 Pack
	1,000	6	<b>60109-1000-6-9</b>	30 Pack
	2,000	12	<b>60109-2000-12-9</b>	20 Pack
	4,000	25	<b>60109-4000-25-9</b>	15 Pack
	20,000	60	<b>60109-20000-60-9</b>	10 Pack

## HyperSep SLE 96 Well Plates

Special treated diatomite SLE	Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
pH=7	200	2	<b>60109-200-2-7W</b>	1 Each
	300	2	<b>60109-300-2-7W</b>	1 Each
	400	2	<b>60109-400-2-7W</b>	1 Each
	500	2	<b>60109-500-2-7W</b>	1 Each
pH=9	200	2	<b>60109-200-2-9W</b>	1 Each
	300	2	<b>60109-300-2-9W</b>	1 Each
	400	2	<b>60109-400-2-9W</b>	1 Each
	500	2	<b>60109-500-2-9W</b>	1 Each

# HyperSep Protein Precipitation Plate

Provides a quick and effective approach for removal of proteins from biological compounds using the protein crash technique. In combination with SPE and SLE, the Protein Precipitation Plate offers a comprehensive range of options for sample preparation of biological based compounds. The 96 well plate format provides a high throughput functionality which lends itself to full automation.

- Dual frit design
- Hydrophobic/oleophobic frits to enable only precipitation of proteins
- Pore size optimized for ideal flow rate
- Specially selected polypropylene for low extractables



## HyperSep Protein Precipitation Plate

Description	Cat. No.	Quantity
Protein Precipitation Plate	60304-201	1 Each



# Microscale SPE Plates and Tips

## HyperSep Lab Plates

For the purification and sample preparation of proteins, DNA, RNA and other biomolecules

- Sample concentration of small-scale samples
- Available in a range of chromatographic materials
- 96 well plate format with media embedded at the bottom of the plate
- Can be processed manually or using a liquid-handling robot
- Not suitable for use with a vacuum

### Applications:

- Tissue culture and separation of products
- Sample concentration
- Sample clean-up
- Collection of sample after chromatography

### HyperSep Lab Plates

Description	Cat. No. (Polystyrene)	Cat. No. (Polypropylene)	Quantity
C18	60110-201	60110-301	5 Pack
C8	60110-202	60110-302	5 Pack
C4	60110-203	60110-303	5 Pack
Hypercarb	60110-204	60110-304	5 Pack
Hypercarb and C18 (Mixed Mode)	60110-205	60110-305	5 Pack
Zirconium Dioxide	60110-206	60110-306	5 Pack
Titanium Dioxide	60110-207	60110-307	5 Pack
SCX	60110-208	60110-308	5 Pack
SAX	60110-209	60110-309	5 Pack

## HyperSep Filter Plates

For effective clean-up of small-scale samples

- 96 well plate for the purification and separation of proteins, peptides, DNA, RNA and other biomolecules
- Clean up of microgram-level samples
- Available in a range of chromatographic materials
- Can be used under vacuum

### HyperSep Filter Plates

Description	Cat. No. (5-7 $\mu$ L Bed Volume)	Cat. No. (40 $\mu$ L Bed Volume)	Quantity
C18	60110-401	60110-501	1 Each
C8	60110-402	60110-502	1 Each
C4	60110-403	60110-503	1 Each
Hypercarb	60110-404	60110-504	1 Each
Hypercarb and C18 (Mix Mode)	60110-405	60110-505	1 Each
Zirconium Dioxide	60110-406	60110-506	1 Each
Titanium Dioxide	60110-407	60110-507	1 Each
SCX	60110-408	60110-508	1 Each
SAX	60110-409	60110-509	1 Each



## HyperSep Microscale Solid Phase Extraction Tips

Revolutionary micropipette tip for sample preparation

- Faster sample preparation with minimal sample loss
- Patented micropipette tip in which the chromatographic material is directly attached to its inner surface
- No contamination from the supporting matrix
- Separation in volumes as low as 100nL
- Applications include mass spectrometry, desalting, protein purification and electrophoresis



### HyperSep Tip Microscale Solid Phase Extraction Tips

Description	Cat. No. (1-10µL Capacity)	Cat. No. (10-200µL Capacity)	Quantity
Thermo Scientific™ BioBasic™ 18	60109-201	60109-209	96 Pack
BioBasic 8	60109-202	60109-210	96 Pack
BioBasic 4	60109-203	60109-211	96 Pack
Hypercarb	60109-204	60109-212	96 Pack
Hypercarb and C18 (mix mode)	60109-205	60109-213	96 Pack
HILIC	60109-206	60109-214	96 Pack
Trypsin	60109-207	60109-215	96 Pack
Titanium Dioxide	60109-208	60109-216	96 Pack
Zirconium Dioxide	60109-217	60109-218	96 Pack

## HyperSep SpinTip Microscale Solid Phase Extraction Tips

Revolutionary micropipette tip for sample preparation

- Pipette tips with a 1 to 2µm wide slit at the bottom that permits the liquid to pass through but retains the chromatographic material (20 to 30µm)
- Faster sample preparation with minimal sample loss
- No contamination from the supporting matrix
- Separation in volumes as low as 100nL



### HyperSep SpinTip Microscale Solid Phase Extraction Tips

Description	Cat. No. (1-10µL Capacity)	Cat. No. (10-200µL Capacity)	Quantity
C18	60109-401	60109-412	96 Pack
C8	60109-402	60109-413	96 Pack
C4	60109-403	60109-414	96 Pack
Hypercarb	60109-404	60109-415	96 Pack
Hypercarb and C18 (mix mode)	60109-405	60109-416	96 Pack
HILIC	60109-406	60109-417	96 Pack
Trypsin	60109-407	60109-418	96 Pack
POROS Weak Anion Exchanger	60109-408	60109-419	96 Pack
POROS Strong Anion Exchanger	60109-409	60109-420	96 Pack
POROS Strong Cation Exchanger	60109-410	60109-421	96 Pack
Titanium Dioxide	60109-411	60109-422	96 Pack
Zirconium Dioxide	60109-424	60109-425	96 Pack

# HyperSep SPE Cartridges and Plates

## HyperSep Retain SPE Cartridges and Plates

Versatile polymeric materials for retention of polar, non-polar, basic and acidic drugs

- Exceptional recoveries for polar and non-polar analytes
- High and consistent recoveries
- High capacity, high purity, highly porous polystyrene DVB material modified with functional groups
  - PEP - urea for polar and non-polar analytes
  - CX - Sulfonic acid for basic analytes
  - AX - Quaternary amine for acidic analytes
- Fast and easy sample preparation and method development
- pH stable 0 to 14



### HyperSep Retain PEP/CX/AX SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Retain PEP Cat No.	Retain CX Cat No.	Retain AX Cat No.	Quantity
30	1	60107-201	60107-301	60107-401	100 Pack
30	3	60107-202	60107-302	60107-402	50 Pack
60	3	60107-203	60107-303	60107-403	50 Pack
60	6	60107-208	60107-308	60107-408	30 Pack
100	3	60107-217	60107-317	60107-417	50 Pack
100	6	60107-207	60107-307	60107-407	30 Pack
150	3	60107-209	60107-309	60107-413	30 Pack
150	6	60107-211	60107-311	60107-411	30 Pack
200	3	60107-204	60107-304	60107-404	50 Pack
200	6	60107-212	60107-314	60107-412	30 Pack
500	3	60107-205	60107-305	60107-405	50 Pack
500	6	60107-206	60107-306	60107-406	30 Pack
500	15	60107-219	60107-319	60107-419	20 Pack
1,000	25	60107-215	60107-315	60107-415	20 Pack
1,000	6	60107-218	60107-318	60107-418	30 Pack
2,000	25	60107-214	60107-312	60107-414	20 Pack

**Note:** For bed weights less than 30mg please refer to SOLA SPE products on page 1-011.

### HyperSep Retain PEP/CX/AX 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Retain PEP Cat No.	Retain CX Cat No.	Retain AX Cat No.	Quantity
30	1	60306-207	60306-303	60306-403	1 Each
60	1	60306-208	60306-304	60306-404	1 Each

**Note:** For bed weights less than 30mg please refer to SOLA SPE products on page 1-011.

### HyperSep Retain PEP/CX/AX 96 Removable Well Plate

Bed Weight (mg)	Volume (mL)	Retain PEP Cat No.	Retain CX Cat No.	Retain AX Cat No.	Quantity
30	1	60303-207	60303-307	60303-407	1 Each
60	1	60303-208	60303-308	60303-408	1 Each

**Note:** For bed weights less than 30mg please refer to SOLA SPE products on page 1-011.

### Retain PEP Applications

- Drugs and metabolites in biological fluids
- Peptides in serum, plasma or biological fluids
- Environmental samples

### Retain CX Applications

- Analysis of a wide range of drugs of abuse including basic and neutral drugs

### Retain AX Applications

- Analysis of THC and its metabolites

## HyperSep Hypercarb SPE Cartridges and Plates

Unique material for retention of highly polar compounds

- Thermo Scientific™ HyperSep™ Hypercarb™ SPE contains flat, 100% porous graphitic carbon (PGC) with selectivity for structurally similar compounds, offering separation of compounds with simple solvents
- pH stable 0 to 14
- High batch-to-batch reproducibility
- Strong retention properties allow use of low bed weights for concentrated extracts
- Interaction mechanism with polar molecules



### HyperSep Hypercarb SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
25	1	<b>60106-304</b>	50 Pack
50	1	<b>60106-303</b>	50 Pack
100	1	<b>60106-302</b>	30 Pack
200	3	<b>60106-301</b>	30 Pack
500	6	<b>60106-402</b>	20 Pack
1,000	6	<b>60106-403</b>	10 Pack
2,000	15	<b>60106-404</b>	10 Pack

### HyperSep Hypercarb 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	<b>60306-501</b>	1 Each
25	1	<b>60306-502</b>	1 Each
50	1	<b>60306-503</b>	1 Each

### HyperSep Hypercarb 96 Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	<b>60302-606</b>	1 Each
25	1	<b>60302-607</b>	1 Each
50	1	<b>60302-608</b>	1 Each

## HyperSep C18 SPE Cartridges and Plates

Feature a highly retentive alkyl-bonded phase for extraction of non-polar to moderately polar compounds

### HyperSep C18 SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
25	1	<b>60108-376</b>	100 Pack
50	1	<b>60108-390</b>	100 Pack
100	1	<b>60108-302</b>	100 Pack
100	3	<b>60108-765</b>	50 Pack
200	3	<b>60108-303</b>	50 Pack
500	3	<b>60108-304</b>	50 Pack
500	6	<b>60108-305</b>	30 Pack
500	10	<b>60108-786</b>	50 Pack
1,000	6	<b>60108-301</b>	30 Pack
1,000	15	<b>60108-776</b>	20 Pack
2,000	15	<b>60108-701</b>	20 Pack
2000	25	<b>60108-780</b>	20 Pack
5,000	25	<b>60108-702</b>	20 Pack
10,000	75	<b>60108-703</b>	10 Pack

### HyperSep C18 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-201</b>	<b>60300-425</b>	1 Each
25	1	<b>60307-202</b>	<b>60300-426</b>	1 Each
50	1	<b>60307-203</b>	<b>60300-427</b>	1 Each
100	1	<b>60307-204</b>	<b>60300-428</b>	1 Each

## HyperSep C8 SPE Cartridges and Plates

Less retentive alternative to C18 for polar and non-polar compounds

### HyperSep C8 SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	<b>60108-391</b>	100 Pack
100	1	<b>60108-392</b>	100 Pack
200	3	<b>60108-393</b>	50 Pack
500	3	<b>60108-309</b>	50 Pack
500	6	<b>60108-394</b>	30 Pack
1,000	6	<b>60108-427</b>	30 Pack
2,000	15	<b>60108-704</b>	20 Pack
5,000	25	<b>60108-705</b>	20 Pack
10,000	75	<b>60108-706</b>	10 Pack

### HyperSep C8 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-211</b>	<b>60300-445</b>	1 Each
25	1	<b>60307-212</b>	<b>60300-446</b>	1 Each
50	1	<b>60307-213</b>	<b>60300-447</b>	1 Each
100	1	<b>60307-214</b>	<b>60300-448</b>	1 Each

### Applications:

- Retentive for non-polar compounds
- Retains most organic analytes from aqueous matrices

### Applications:

- Drugs and their metabolites in biological samples
- Peptides in biological samples

## HyperSep Silica SPE Cartridges and Plates

A polar sorbent primarily used to retain analytes in non-polar matrices

### HyperSep Silica SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	<b>60108-409</b>	100 Pack
100	1	<b>60108-317</b>	100 Pack
200	3	<b>60108-410</b>	50 Pack
500	3	<b>60108-315</b>	50 Pack
500	6	<b>60108-411</b>	30 Pack
500	10	<b>60108-793</b>	50 Pack
1,000	6	<b>60108-426</b>	30 Pack
2,000	15	<b>60108-710</b>	20 Pack
5,000	25	<b>60108-711</b>	20 Pack
10,000	25	<b>60108-853</b>	20 Pack
10,000	75	<b>60108-712</b>	10 Pack
50,000	150	<b>60108-850</b>	10 Pack
20,000	75	<b>60108-851</b>	10 Pack
70,000	150	<b>60108-852</b>	10 Pack

### HyperSep Silica 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-231</b>	<b>60300-485</b>	1 Each
25	1	<b>60307-232</b>	<b>60300-486</b>	1 Each
50	1	<b>60307-233</b>	<b>60300-487</b>	1 Each
100	1	<b>60307-234</b>	<b>60300-488</b>	1 Each

### Applications:

- Extraction of polar compounds including aldehydes, amines, drugs, pesticides and herbicides
- Extraction of carotenoids, fat-soluble vitamins, aflatoxins in food matrices
- Extraction of fatty acids and phospholipids

## HyperSep Acid Wash Silica SPE Cartridges and Plates

### HyperSep Acid Wash Silica SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	<b>60108-901</b>	100 Pack
100	1	<b>60108-902</b>	100 Pack
100	10	<b>60108-903</b>	50 Pack
200	3	<b>60108-904</b>	50 Pack
200	10	<b>60108-905</b>	50 Pack
500	3	<b>60108-906</b>	50 Pack
500	6	<b>60108-907</b>	30 Pack
500	10	<b>60108-908</b>	50 Pack
1,000	6	<b>60108-909</b>	30 Pack
2,000	15	<b>60108-910</b>	20 Pack
5,000	25	<b>60108-911</b>	20 Pack
10,000	75	<b>60108-912</b>	10 Pack

### Applications:

- Extraction of polar compounds from ground, waste and drinking water
- Extraction of polar compounds from food and beverages

### HyperSep Acid Wash Silica 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-321</b>	<b>60300-831</b>	1 Each
25	1	<b>60307-322</b>	<b>60300-832</b>	1 Each
50	1	<b>60307-323</b>	<b>60300-833</b>	1 Each
100	1	<b>60307-324</b>	<b>60300-834</b>	1 Each

## HyperSep Phenyl SPE Cartridges and Plates

Offer alternative selectivity for retention of basic compounds

### HyperSep Phenyl SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	<b>60108-516</b>	100 Pack
100	1	<b>60108-386</b>	100 Pack
200	3	<b>60108-387</b>	50 Pack
500	3	<b>60108-388</b>	50 Pack
500	6	<b>60108-389</b>	30 Pack
1,000	6	<b>60108-517</b>	30 Pack
2,000	15	<b>60108-707</b>	20 Pack
5,000	25	<b>60108-708</b>	20 Pack
10,000	75	<b>60108-709</b>	10 Pack

### Applications:

- Extraction of aromatic compounds
- Extraction of basic analytes

### HyperSep Phenyl 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-221</b>	<b>60300-685</b>	1 Each
25	1	<b>60307-222</b>	<b>60300-686</b>	1 Each
50	1	<b>60307-223</b>	<b>60300-687</b>	1 Each
100	1	<b>60307-224</b>	<b>60300-688</b>	1 Each

## HyperSep Diol SPE Cartridges and Plates

Ideal for the extraction of polar compounds

### HyperSep Diol SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	<b>60108-571</b>	100 Pack
100	1	<b>60108-572</b>	100 Pack
200	3	<b>60108-573</b>	50 Pack
500	3	<b>60108-574</b>	50 Pack
500	6	<b>60108-575</b>	30 Pack
1,000	6	<b>60108-576</b>	30 Pack
2,000	15	<b>60108-755</b>	20 Pack
5,000	25	<b>60108-756</b>	20 Pack
10,000	75	<b>60108-757</b>	10 Pack

### Applications:

- Normal phase extraction of polar compounds
- Purification of polar compounds

### HyperSep Diol 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-311</b>	<b>60300-630</b>	1 Each
25	1	<b>60307-312</b>	<b>60300-631</b>	1 Each
50	1	<b>60307-313</b>	<b>60300-632</b>	1 Each
100	1	<b>60307-314</b>	<b>60300-633</b>	1 Each

## HyperSep Verify CX SPE Cartridges and Plates

Features non-polar and ionic characteristics for improved analysis of drugs of abuse

### HyperSep Verify CX SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	<b>60108-741</b>	100 Pack
130	1	<b>60108-719</b>	100 Pack
130	10	<b>60108-769</b>	50 Pack
200	3	<b>60108-777</b>	50 Pack
200	6	<b>60108-722</b>	50 Pack
200	10	<b>60108-742</b>	50 Pack
300	3	<b>60108-720</b>	50 Pack
500	3	<b>60108-721</b>	50 Pack
500	6	<b>60108-723</b>	30 Pack
1,000	6	<b>60108-724</b>	30 Pack

#### Applications:

- Analysis of a wide range of drugs of abuse from biological matrices, including basic and neutral drugs

### HyperSep Verify CX 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-261</b>	<b>60300-805</b>	1 Each
25	1	<b>60307-262</b>	<b>60300-806</b>	1 Each
50	1	<b>60307-263</b>	<b>60300-807</b>	1 Each
100	1	<b>60307-264</b>	<b>60300-808</b>	1 Each

## HyperSep Verify AX SPE Cartridges and Plates

Features non-polar and ionic characteristics for improved analysis of acidic drugs and metabolites

### HyperSep Verify AX SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
130	1	<b>60108-727</b>	100 Pack
130	10	<b>60108-767</b>	50 Pack
200	3	<b>60108-768</b>	30 Pack
200	6	<b>60108-730</b>	50 Pack
200	10	<b>60108-764</b>	50 Pack
300	3	<b>60108-728</b>	50 Pack
500	3	<b>60108-729</b>	50 Pack
500	6	<b>60108-731</b>	30 Pack
1,000	6	<b>60108-732</b>	30 Pack

#### Applications:

- Analysis of THC and its metabolites

### HyperSep Verify AX 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-271</b>	<b>60300-813</b>	1 Each
25	1	<b>60307-272</b>	<b>60300-814</b>	1 Each
50	1	<b>60307-273</b>	<b>60300-815</b>	1 Each
100	1	<b>60307-274</b>	<b>60300-816</b>	1 Each

## HyperSep Strong Cation Exchanger (SCX) SPE Cartridges and Plates

A strong cation exchange sorbent for extraction of charged basic compounds

### HyperSep Strong Cation Exchanger (SCX) SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	<b>60108-420</b>	100 Pack
100	1	<b>60108-421</b>	100 Pack
200	3	<b>60108-422</b>	50 Pack
500	3	<b>60108-423</b>	50 Pack
500	6	<b>60108-520</b>	30 Pack
1,000	6	<b>60108-433</b>	30 Pack
2,000	15	<b>60108-716</b>	20 Pack
5,000	25	<b>60108-717</b>	20 Pack
10,000	75	<b>60108-718</b>	10 Pack
20,000	75	<b>60108-857</b>	10 Pack
50,000	150	<b>60108-858</b>	10 Pack

#### Applications:

- Isolation of cationic compounds from complex samples
- Removal of basic drugs
- Organic bases
- Nucleic acid bases, nucleotides and surfactants

### HyperSep Strong Cation Exchanger (SCX) 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-251</b>	<b>60300-585</b>	1 Each
25	1	<b>60307-252</b>	<b>60300-586</b>	1 Each
50	1	<b>60307-253</b>	<b>60300-587</b>	1 Each
100	1	<b>60307-254</b>	<b>60300-588</b>	1 Each

## HyperSep Strong Anion Exchanger (SAX) SPE Cartridges and Plates

Strong anion exchange sorbent for extraction of weak acids

### HyperSep Strong Anion Exchanger (SAX) SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	<b>60108-417</b>	100 Pack
100	1	<b>60108-418</b>	100 Pack
200	3	<b>60108-419</b>	50 Pack
500	3	<b>60108-521</b>	50 Pack
500	6	<b>60108-360</b>	30 Pack
1,000	6	<b>60108-434</b>	30 Pack
2,000	15	<b>60108-713</b>	20 Pack
5,000	25	<b>60108-714</b>	20 Pack
10,000	75	<b>60108-715</b>	10 Pack

#### Applications:

- Isolation of anionic compounds from complex samples
- Removal of acidic food pigments
- Isolation of phenolic compounds
- Nucleic acids, nucleotides and surfactants

### HyperSep Strong Anion Exchanger (SAX) 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-241</b>	<b>60300-565</b>	1 Each
25	1	<b>60307-242</b>	<b>60300-566</b>	1 Each
50	1	<b>60307-243</b>	<b>60300-567</b>	1 Each
100	1	<b>60307-244</b>	<b>60300-568</b>	1 Each



## HyperSep Carboxylic Acid (WCX) SPE Cartridges and Plates

Weak cation exchange sorbent for extraction of strong bases

### HyperSep Carboxylic Acid (WCX) SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	<b>60108-201</b>	100 Pack
100	1	<b>60108-202</b>	100 Pack
200	3	<b>60108-203</b>	50 Pack
200	10	<b>60108-207</b>	50 Pack
500	3	<b>60108-204</b>	50 Pack
500	6	<b>60108-205</b>	30 Pack
1,000	6	<b>60108-206</b>	30 Pack
2,000	15	<b>60108-208</b>	20 Pack
5,000	25	<b>60108-209</b>	20 Pack

#### Applications:

- Drugs and metabolites in biological fluids
- Peptides in serum, plasma or biological fluids
- Pesticides in water, food and beverages
- Nucleic acids, nucleotides and surfactants

### HyperSep Carboxylic Acid (WCX) 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-331</b>	<b>60300-841</b>	1 Each
25	1	<b>60307-332</b>	<b>60300-842</b>	1 Each
50	1	<b>60307-333</b>	<b>60300-843</b>	1 Each
100	1	<b>60307-334</b>	<b>60300-844</b>	1 Each

## HyperSep Aminopropyl (WAX) SPE Cartridges and Plates

Weak anion exchange sorbent for extraction of strong acids

### HyperSep Aminopropyl (WAX) SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	<b>60108-424</b>	100 Pack
50	3	<b>60108-429</b>	50 Pack
100	1	<b>60108-364</b>	100 Pack
200	3	<b>60108-425</b>	50 Pack
500	3	<b>60108-518</b>	50 Pack
500	6	<b>60108-519</b>	30 Pack
1,000	6	<b>60108-432</b>	30 Pack
2,000	15	<b>60108-738</b>	20 Pack
5,000	25	<b>60108-739</b>	20 Pack
10,000	75	<b>60108-740</b>	10 Pack

#### Applications:

- Separation of structural isomers
- Drugs and metabolites in biological fluids
- Separation of saccharides, phenols and petroleum products

### HyperSep Aminopropyl (WAX) 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-291</b>	<b>60300-505</b>	1 Each
25	1	<b>60307-292</b>	<b>60300-506</b>	1 Each
50	1	<b>60307-293</b>	<b>60300-507</b>	1 Each
100	1	<b>60307-294</b>	<b>60300-508</b>	1 Each

## HyperSep Florisil SPE Cartridges and Plates

Ideal for the isolation of polar compounds from non-polar matrices

### HyperSep Florisil SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	<b>60108-402</b>	100 Pack
100	1	<b>60108-403</b>	100 Pack
200	3	<b>60108-404</b>	50 Pack
500	3	<b>60108-405</b>	50 Pack
500	6	<b>60108-500</b>	30 Pack
1,000	6	<b>60108-431</b>	30 Pack
2,000	6	<b>60108-758</b>	30 Pack
2,000	15	<b>60108-735</b>	20 Pack
5,000	25	<b>60108-736</b>	20 Pack
10,000	75	<b>60108-737</b>	10 Pack

### HyperSep Florisil 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-281</b>	<b>60300-725</b>	1 Each
25	1	<b>60307-282</b>	<b>60300-726</b>	1 Each
50	1	<b>60307-283</b>	<b>60300-727</b>	1 Each
100	1	<b>60307-284</b>	<b>60300-728</b>	1 Each

## HyperSep Cyano SPE Cartridges and Plates

Optimized for the retention of moderately polar compounds from non-polar matrices

### HyperSep Cyano SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	<b>60108-746</b>	100 Pack
100	1	<b>60108-745</b>	100 Pack
200	3	<b>60108-747</b>	50 Pack
500	3	<b>60108-748</b>	50 Pack
500	6	<b>60108-749</b>	30 Pack
1,000	6	<b>60108-750</b>	30 Pack
2,000	15	<b>60108-751</b>	20 Pack
5,000	25	<b>60108-752</b>	20 Pack
10,000	75	<b>60108-753</b>	10 Pack

### HyperSep Cyano 96 Fixed/Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Fixed Cat. No.	Removable Cat. No.	Quantity
10	1	<b>60307-301</b>	<b>60300-821</b>	1 Each
25	1	<b>60307-302</b>	<b>60300-822</b>	1 Each
50	1	<b>60307-303</b>	<b>60300-823</b>	1 Each
100	1	<b>60307-304</b>	<b>60300-824</b>	1 Each



### Applications:

- Extraction of pesticides using official methods
- Polychlorinated biphenyls in transformer oil
- Alcohol, aldehydes, amines and drugs

### Applications:

- Retaining polar compounds from hexane and oils
- Reversed phase extraction of moderately polar compounds

## SPE Manifolds and Accessories

A comprehensive range of accessories to improve your sample preparation workflow

### HyperSep Positive Pressure Manifold

#### HyperSep Positive Pressure Manifold

Description	Cat. No.	Quantity
Positive Pressure Manifold with 13mm Collection Rack	<b>60104-236</b>	1 Each
Positive Pressure Manifold with 16mm Collection Rack	<b>60104-274</b>	1 Each
Adaptor Plate for 1mL Cartridges for PP Manifold	<b>60104-265</b>	1 Each
Adaptor Plate for 3mL Cartridges for PP Manifold	<b>60104-266</b>	1 Each
Adaptor Plate for 6mL Cartridges for PP Manifold	<b>60104-267</b>	1 Each
Adaptor Plate for 10mL/15mL Cartridges for PP Manifold	<b>60104-271</b>	1 Each
Collection Rack for 13mm Tubes	<b>60104-268</b>	1 Each
Collection Rack for 16mm Tubes	<b>60104-269</b>	1 Each
Pre-Drilled Waste Container	<b>60104-270</b>	1 Each
Installation Kit	<b>60104-272</b>	1 Each
In Line Air Filter	<b>60104-273</b>	1 Each



## HyperSep Universal Vacuum Manifold

- Accommodates both SPE cartridges and 96 well plates
- System supplied with: Manifold, base/gauge, flask and stopper, tubing and spigots




### HyperSep Universal Vacuum Manifold

Description	Cat. No.	Quantity
Universal Vacuum Manifold	<b>60104-230</b>	1 Each
Vacuum Pump, European Version	<b>60104-241</b>	1 Each
Vacuum Pump, North American Version	<b>60104-243</b>	1 Each
24 Well Extraction Plate	<b>60104-284</b>	1 Each
48 Well Extraction Plate	<b>60104-289</b>	1 Each
Plugs for 24-position Extraction Plate	<b>60104-234</b>	24 Pack
Plugs for 48-position Extraction Plate	<b>60104-235</b>	48 Pack

## HyperSep-96 Well Plate Vacuum Manifold

Included with system:

- Base
- Lid
- Waste collection tray



HyperSep-96 well plate is not included and needs to be purchased separately.

### HyperSep-96 Well Plate Manifolds

Description	Cat. No.	Quantity
HyperSep-96 Vacuum Manifold	<b>60103-351</b>	1 Each
Vacuum Pump, European Version	<b>60104-241</b>	1 Each
Vacuum Pump, North American Version	<b>60104-243</b>	1 Each


## HyperSep Glass Block Manifolds

16-Port Vacuum Manifold:

- Glass Block, Corian Manifold Lid, Cover Gasket, Vacuum Gauge and Valve Assembly, 16 Tips, Adjustable Collection Rack, Bulkhead Luer Fittings, 16 Plugs and Manifold Safety Tray

24-Port Vacuum Manifold:

- Glass Block, Corian Manifold Lid, Cover Gasket, Vacuum Gauge and Valve Assembly, 24 Tips, Adjustable Collection Rack, Bulkhead Luer Fittings, 24 Plugs and Manifold Safety Tray



Luer Lock Plugs  
(Cat. No. 60104-258)



Retaining Clips for Collection Rack  
(Cat. No. 60104-255)



Bulkhead Luer Fittings  
(Cat. No. 60104-256)



Collection Rack  
(Cat. No. 60104-251 (16 port),  
60104-252 (24 port))



Vacuum Gauge  
(Cat. No. 60104-240)



Valve Assembly  
(Cat. No. 60104-261)




Tips  
(Cat. No. 60104-245)



Manifold Lid Legs  
(Cat. No. 60104-257)

## HyperSep Glass Block Manifolds

Description	Cat. No.	Quantity
16-port Vacuum Manifold	<b>60104-232</b>	1 Each
24-port Vacuum Manifold	<b>60104-233</b>	1 Each
Vacuum Pump, European Version	<b>60104-241</b>	1 Each
Vacuum Pump, North American Version	<b>60104-243</b>	1 Each
<b>Replacement Parts</b>		
Vacuum Gauge	<b>60104-240</b>	1 Each
Stopcocks for 16-port Vacuum Manifold	<b>60104-242</b>	16 Pack
Stopcocks for 24-port Vacuum Manifold	<b>60104-244</b>	24 Pack
Tips for Vacuum Manifold	<b>60104-245</b>	12 Pack
Lid for 16-port Glass Block Manifold	<b>60104-262</b>	1 Each
Lid for 24-port Glass Block Manifold	<b>60104-248</b>	1 Each
Gasket for 16-port Manifold	<b>60104-249</b>	1 Each
Gasket for 24-port Manifold	<b>60104-250</b>	1 Each
Collection Rack for 16-port Vacuum Manifold	<b>60104-251</b>	1 Each
Collection Rack for 24-port Vacuum Manifold	<b>60104-252</b>	1 Each
Glass Block for 16-port Vacuum Manifold	<b>60104-253</b>	1 Each
Glass Block for 24-port Vacuum Manifold	<b>60104-254</b>	1 Each
Retaining Clips for Collection Rack	<b>60104-255</b>	12 Pack
Bulkhead Luer Fittings	<b>60104-256</b>	12 Pack
Manifold Lid Legs	<b>60104-257</b>	4 Pack
Luer Lock Plugs	<b>60104-258</b>	12 Pack
Column Adaptors	<b>60104-259</b>	15 Pack
Manifold Safety Tray	<b>60104-260</b>	1 Each
Vacuum Gauge and Valve Assembly	<b>60104-261</b>	1 Each
Collection Rack Posts	<b>60104-264</b>	3 Pack



## HyperSep Online SPE Products

Retain specific analytes in a sample matrix when used with an appropriate HPLC column



- Effective removal of contaminants such as proteins from samples
- Compatible with conventional HPLC systems
- Fast and effective clean-up and concentration of target compounds

### HyperSep Javelin Direct-Connect Online SPE Columns

ID (mm)	Length (mm)	Retain PEP	Retain CX	Retain AX	Hypercarb	Quantity
2.1	10	<b>60310-201</b>	<b>60310-301</b>	<b>60310-401</b>	<b>60310-501</b>	4 Pack
3.0	10	<b>60310-202</b>	<b>60310-302</b>	<b>60310-402</b>	<b>60310-502</b>	4 Pack

### HyperSep UNIGUARD Direct-Connect Online SPE Cartridges

ID (mm)	Length (mm)	Retain PEP	Retain CX	Retain AX	Hypercarb	Quantity
2.1	10	<b>60311-201</b>	<b>60311-301</b>	<b>60311-401</b>	<b>60311-501</b>	4 Pack
3.0	10	<b>60311-202</b>	<b>60311-302</b>	<b>60311-402</b>	<b>60311-502</b>	4 Pack

### HyperSep HPLC Columns for Online SPE

ID (mm)	Length (mm)	Retain PEP	Retain CX	Retain AX	Hypercarb	Quantity
2.1	20	<b>60312-201</b>	<b>60312-301</b>	<b>60312-401</b>	<b>60312-501</b>	1 Each
3.0	20	<b>60312-202</b>	<b>60312-302</b>	<b>60312-402</b>	<b>60312-502</b>	1 Each

*Perform Better. Live Longer.*

*Discover our Accucore Solid Core LC Columns*

» Visit **PAGE 4-005**





## HyperSep Dispersive SPE (QuEChERS)

### Fast, easy and robust sample preparation of food for pesticide analysis

QuEChERS dispersive SPE products provide a Quick, Easy, Cheap, Effective, Rugged and Safe sample preparation technique for multi-residue pesticide analysis. Products are available in a wide range of formats to meet application requirements.

QuEChERS methods overcome problems associated with time consuming, expensive and labor intensive multi-residue methods (MRM).

- **Fast and Easy Sample Preparation**
- **Robust and Reproducible Process**
- **High Recovery and Reproducibility**



Learn more about HyperSep Dispersive SPE at  
[www.thermoscientific.com/samplepreparation](http://www.thermoscientific.com/samplepreparation)



## HyperSep Dispersive SPE

### Considerations in method development

- Determine the properties of the pesticides of interest:
  - Base sensitive compounds - use extraction product with sodium acetate
  - pH dependent
  - Non-base sensitive compounds - use extraction product with sodium chloride or sodium citrate
- Determine the properties of the sample matrix:
  - General
  - Fatty
  - Pigmented
  - Highly pigmented

### Product Selection – Sample Cleanup

Matrix Type	Examples	Sorbent Requirements
General Matrices	Apples Cucumber Melon	Magnesium sulfate, PSA
Fatty Matrices	Milk Cereals Fish	Magnesium sulfate, PSA, C18
Pigmented Matrices	Lettuce Carrot Wine	Magnesium sulfate, PSA, C18, GCB
High Pigmented Matrices	Spinach Red Peppers	Magnesium sulfate, PSA, C18, GCB



### QuEChERS Troubleshooting

Problem	Causes	Recommended Solutions
Loss of planar pesticides	Presence of GCB may result in a loss of planar compounds	Use a product with less GCB Use the Dual Phase QuEChERS product
Loss of acidic compounds e.g. 2,4-D from starting matrix	Presence of PSA will extract acidic compounds from matrix	Use a product containing magnesium sulfate and C18
Loss of compounds during subsequent analysis	Some compounds are unstable and can break down during analysis	Use an analyte protectant e.g. toluene or sorbitol
Addition of sample to QuEChERS extraction tube containing sorbent causes an exothermic reaction	Exothermic reaction between water in sample and magnesium sulfate	Add the sample to the tube, then the solvents, then the sorbent materials
Poor recovery of pesticide compounds	Sample not in appropriate homogenization state	Wrong products used in method Ensure sample is hydrated to 80% or higher Verify nature of pesticides e.g. are base sensitive compounds present

# HyperSep Dispersive SPE Extraction Products

## HyperSep Dispersive SPE Pre-Packed Extraction Products

Description	Capacity (mL)	Cat. No. (250 Pack)	Cat. No. (25 Pack)
6000mg anhydrous magnesium sulfate, 1500mg anhydrous sodium acetate	50	<b>60105-210</b>	<b>60105-310</b>
4000mg anhydrous magnesium sulfate, 1000mg anhydrous sodium chloride	50	<b>60105-211</b>	<b>60105-311</b>
4000mg anhydrous magnesium sulfate, 1000mg sodium chloride, 1000mg sodium citrate tribasic dihydrate, 500mg disodium citrate sesquihydrate	50	<b>60105-216</b>	<b>60105-316</b>
6000mg anhydrous magnesium sulfate, 1500mg sodium chloride, 1500mg sodium citrate tribasic dihydrate, 750mg disodium citrate sesquihydrate	50	<b>60105-212</b>	<b>60105-312</b>
1500mg magnesium sulfate, 500mg C18	50	<b>60105-234</b>	–

**Each pack contains centrifuge tubes pre-packed with sorbents**

## HyperSep Dispersive SPE Multipacks

Description	Cat. No.	Quantity
4000mg anhydrous magnesium sulfate, 1000mg sodium chloride	<b>60105-332</b>	50 Pack
4000mg anhydrous magnesium sulfate, 1000mg sodium chloride, 500mg sodium citrate dibasic sesquihydrate, 1000mg sodium citrate tribasic	<b>60105-333</b>	50 Pack
4000mg anhydrous magnesium sulfate, 1000mg sodium acetate	<b>60105-334</b>	50 Pack
6000mg anhydrous magnesium sulfate, 1500mg sodium acetate	<b>60105-335</b>	50 Pack
6000mg anhydrous magnesium sulfate, 1500mg sodium chloride	<b>60105-336</b>	50 Pack
6000mg anhydrous magnesium sulfate, 1500mg magnesium sulfate, 1500mg sodium citrate dihydrate, 750mg disodium citrate sesquihydrate	<b>60105-337</b>	50 Pack
8000mg anhydrous magnesium sulfate, 2000mg of sodium chloride	<b>60105-338</b>	50 Pack
8000mg anhydrous magnesium sulfate, 3500mg of sodium chloride	<b>60105-339</b>	50 Pack

**Each pack contains centrifuge tubes and sorbents packaged in mylar pouches**

## HyperSep Dispersive SPE Mylar Pouches

Description	Cat. No.	Quantity
4000mg anhydrous magnesium sulfate, 1000mg sodium chloride	<b>60105-340</b>	50 Pack
6000mg anhydrous magnesium sulfate, 1500mg sodium acetate (anhydrous)	<b>60105-341</b>	50 Pack
6000mg anhydrous magnesium sulfate, 1500mg sodium chloride	<b>60105-342</b>	50 Pack
6000mg anhydrous magnesium sulfate, 1500mg sodium chloride, 1500mg sodium citrate dihydrate, 750mg disodium citrate sesquihydrate	<b>60105-343</b>	50 Pack
4000mg anhydrous magnesium sulfate, 1000mg sodium chloride, 500mg sodium citrate dibasic sesquihydrate, 1000mg sodium citrate tribasic	<b>60105-344</b>	50 Pack
1200mg sodium chloride	<b>60105-345</b>	50 Pack
4000mg sodium sulfate and 500mg anhydrous magnesium sulfate	<b>60105-346</b>	50 Pack
4000mg anhydrous magnesium sulfate, 1000mg sodium chloride, 1500mg sodium citrate dibasic dihydrate	<b>60105-347</b>	100 Pack
4000mg anhydrous magnesium sulfate, 2000mg sodium chloride, 1500mg sodium citrate dibasic dihydrate	<b>60105-348</b>	100 Pack

**Each pack contains sorbents packaged in mylar pouches (no centrifuge tubes provided)**

## HyperSep Dispersive SPE Clean-Up Products

### HyperSep Dispersive SPE Clean-Up Products

Description	Capacity (mL)	Cat. No. (100 Pack)	Cat. No. (10 Pack)
150mg anhydrous magnesium sulfate, 50mg PSA, 50mg GCB	2	60105-202	60105-302
150mg anhydrous magnesium sulfate, 50mg PSA	2	60105-203	60105-303
150mg anhydrous magnesium sulfate, 50mg PSA, 50mg C18	2	60105-204	60105-304
150mg anhydrous magnesium sulfate, 25mg PSA	2	60105-219	60105-319
150mg anhydrous magnesium sulfate, 25mg PSA, 25mg C18	2	60105-220	60105-320
150mg anhydrous magnesium sulfate, 25mg PSA, 2.5mg GCB	2	60105-221	60105-321
150mg anhydrous magnesium sulfate, 25mg PSA, 7.5mg GCB	2	60105-222	60105-322
150mg anhydrous magnesium sulfate, 50mg PSA, 50mg C18, 50mg GCB	2	60105-223	60105-323
150mg anhydrous magnesium sulfate, 50mg PSA, 50mg chlorofiltr™	2	60105-232	-
150mg anhydrous magnesium sulfate, 150mg PSA	2	-	60105-350

Description	Capacity (mL)	Cat. No. (100 Pack)	Cat. No. (10 Pack)
900mg anhydrous magnesium sulfate, 300mg PSA, 150mg GCB	15	60105-205	60105-305
900mg anhydrous magnesium sulfate, 300mg PSA, 150mg C18	15	60105-206	60105-306
750mg anhydrous magnesium sulfate, 250mg PSA, 250mg C18, 250mg GCB	15	60105-213	60105-313
900mg anhydrous magnesium sulfate, 300mg PSA	15	60105-214	60105-314
900mg anhydrous magnesium sulfate, 150mg PSA	15	60105-215	60105-315
900mg anhydrous magnesium sulfate, 150mg PSA, 45mg GCB	15	60105-217	60105-317
900mg anhydrous magnesium sulfate, 150mg PSA, 15mg GCB	15	60105-218	60105-318
1200mg anhydrous magnesium sulfate, 400mg PSA	15	60105-224	60105-324
1200mg anhydrous magnesium sulfate, 400mg PSA, 400mg C18	15	60105-225	60105-325
1200mg anhydrous magnesium sulfate, 400mg PSA, 400mg C18, 400mg GCB	15	60105-226	60105-326
900mg anhydrous magnesium sulfate, 150mg PSA, 150mg C18	15	60105-227	60105-327
150mg anhydrous magnesium sulfate, 50mg PSA, 50mg GCB	15	60105-230	60105-330
150mg anhydrous magnesium sulfate, 300mg PSA, 150mg chlorofiltr	15	60105-231	60105-331
1200mg anhydrous magnesium sulfate, 900mg PSA	15	60105-235	-
4000mg anhydrous magnesium sulfate, 1000mg sodium chloride, 1000mg sodium citrate tribasic dihydrate, 500mg sodium citrate dibasic sesquihydrate	15	60105-236	-
150mg anhydrous magnesium sulfate, 150mg PSA, 50mg C18	15	60105-238	-

**Each pack contains centrifuge tubes pre-packed with sorbents**

### Homogenizers

Description	Cat. No.	Quantity
Homogenizers for 15mL QuEChERS tubes	60106-421	50 Pack
Homogenizers for 50mL QuEChERS tubes	60106-422	50 Pack

### HyperSep Dispersive SPE Clean-Up Products for Schenck Method

Description	Capacity (mL)	Cat. No.	Quantity
400mg PSA on bottom, 200mg GCB on top, separated by a frit	6	60105-207	30 Pack
500mg PSA on bottom, 250mg GCB on top, separated by a frit	6	60105-208	30 Pack
500mg PSA on bottom, 500mg GCB on top, separated by a frit	6	60105-209	30 Pack



## Syringes and Filtration

### Providing high quality filtration solutions for a variety of samples and applications

Choose from a variety of syringe filter sizes, membranes and housings for a wide range of laboratory applications.

- **High Quality Filtration**  
Low extractable membranes and solvent-resistant housing
- **Titan3 Color Coded Ring**  
For easy selection of correct membrane and pore size
- **High Burst Pressures**  
For enhanced robustness in sample processing



The syringe filter web page contains the latest news, applications and downloads for the range.  
Visit [www.thermoscientific.com/syringefilters](http://www.thermoscientific.com/syringefilters)

# Syringes and Filtration

## Technical Resources



visit [www.thermoscientific.com/chromexpert](http://www.thermoscientific.com/chromexpert)  
to access the following information:

- Membrane Selection Guide
- Membrane Compatibility Chart



## Titan3 Syringe Filters

- Low extractable membranes and housing
- HPLC performance tested
- Color coding for easy selection of the correct membrane and pore size
- Enhanced Luer Lock inlet which prevents leakage
- Most 30mm devices are provided with a 1mm boro-silicate glass pre-filter  
This is of benefit for high solids samples with larger size particulates
- Integral ring provides greater strength to the housing preventing leakage and bursting
- 30mm products pressure rated to 120psi
- Packed in re-usable rigid transparent color coded containers



### Titan3 Nylon Syringe Filters

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
4	0.2	No	<b>42204-NN</b>	100 Pack
	0.45	No	<b>44504-NN</b>	100 Pack
17	0.2	No	<b>42213-NN</b>	200 Pack
	0.45	No	<b>44513-NN</b>	200 Pack
30	0.2	Yes	<b>42225-NN</b>	100 Pack
	0.45	Yes	<b>44525-NN</b>	100 Pack
	0.45	No	<b>44526-NN</b>	100 Pack
	1.5	No	<b>41225-NN</b>	100 Pack
	5	No	<b>45025-NN</b>	100 Pack

### Titan3 PTFE (Hydrophobic) Syringe Filters

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
4	0.2	No	<b>42204-NP</b>	100 Pack
	0.45	No	<b>44504-NP</b>	100 Pack
17	0.2	No	<b>42213-NP</b>	200 Pack
	0.45	No	<b>44513-NP</b>	200 Pack
30	0.2	Yes	<b>42225-NP</b>	100 Pack
	0.45	Yes	<b>44525-NP</b>	100 Pack
	1.0	No	<b>41025-NP</b>	100 Pack

### Titan3 PTFE (Hydrophilic) Syringe Filters

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
17	0.2	No	<b>42213-NPL</b>	200 Pack
	0.45	No	<b>44513-NPL</b>	200 Pack
30	0.2	No	<b>42225-NPL</b>	100 Pack
	0.45	No	<b>44525-NPL</b>	100 Pack

### Titan3 PVDF Syringe Filters

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
4	0.2	No	<b>42204-PV</b>	100 Pack
	0.45	No	<b>44504-PV</b>	100 Pack
17	0.2	Yes	<b>42213-PV</b>	200 Pack
	0.45	No	<b>44513-PV</b>	200 Pack
30	0.2	Yes	<b>42225-PV</b>	100 Pack
	0.45	No	<b>44525-PV</b>	100 Pack

#### Applications:

- General laboratory filtration
- Filtration for most HPLC samples

#### Applications:

- Filtration of aggressive organic, highly basic or hot solutions, transducer protectors
- Filter aqueous solutions after prewetting with an alcohol

#### Applications:

- Filtration of aggressive organic, highly basic or hot solutions, without the need pre-wet membrane
- Filter aqueous solutions without prewetting with an alcohol

#### Applications:

- General biological filtration
- Filtration of samples where high protein recovery is desired

**Titan3 Regenerated Cellulose Syringe Filters**

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
4	0.2	No	<b>52204-RC</b>	100 Pack
	0.45	No	<b>54504-RC</b>	100 Pack
17	0.2	Yes	<b>52213-RC</b>	200 Pack
	0.45	No	<b>54513-RC</b>	200 Pack
30	0.2	Yes	<b>52225-RC</b>	100 Pack
	0.45	No	<b>54525-RC</b>	100 Pack

**Titan3 PES (Polyethersulfone) Syringe Filters**

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
17	0.2	No	<b>42213-PS</b>	200 Pack
	0.45	No	<b>44513-PS</b>	200 Pack
30	0.2	Yes	<b>42225-PS</b>	100 Pack
	0.45	No	<b>44525-PS</b>	100 Pack

**Titan3 GMF (Glass Microfiber) Syringe Filters**

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
30	0.7	No	<b>40725-GM</b>	100 Pack
	1.2	No	<b>41225-GM</b>	100 Pack
	3.1	No	<b>42725-GM</b>	100 Pack

**Titan3 Polypropylene Syringe Filters**

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
17	0.2	No	<b>42213-PP</b>	200 Pack
	0.45	No	<b>44513-PP</b>	200 Pack
30	0.2	No	<b>42225-PP</b>	100 Pack
	0.45	No	<b>44525-PP</b>	100 Pack

**Titan3 Cellulose Acetate Syringe Filters**

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
4	0.2	No	<b>42204-CA</b>	100 Pack
	0.45	No	<b>44502-CA</b>	100 Pack
17	0.2	No	<b>42213-CA</b>	100 Pack
	0.45	No	<b>44513-CA</b>	100 Pack
30	0.2	No	<b>42225-CA</b>	100 Pack
	0.45	No	<b>44525-CA</b>	100 Pack

**Applications:**

- Low non-specific binding applications
- Tissue culture media filtration and general biological sample filtration

**Applications:**

- Ion chromatography
- Tissue culture filtration, filtration of proteins and nucleic acids
- High-temperature liquids

**Applications:**

- Dissolution testing
- General filtration

**Applications:**

- Filtration of biological samples
- Filtration of aggressive organic solutions

**Applications:**

- Tissue culture media filtration, sensitive biological samples



## Target2 Syringe Filters

- Low extractable membranes and housing
- HPLC performance tested
- Plain polypropylene housing
- 30mm products pressure rated to 100psi
- Enhanced Luer Lock inlet which prevents leakage



### Target2 Nylon Syringe Filters

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
4	0.2	No	<b>F2504-2</b>	100 Pack
	0.45	No	<b>F2504-1</b>	100 Pack
17	0.2	No	<b>F2513-2</b>	100 Pack
	0.45	No	<b>F2513-1</b>	100 Pack
30	0.2	No	<b>F2500-2</b>	100 Pack
		Yes	<b>F2502-2</b>	100 Pack
	0.45	No	<b>F2500-1</b>	100 Pack
		Yes	<b>F2502-1</b>	100 Pack
	1.5	No	<b>F2500-12</b>	100 Pack
	5.0	No	<b>F2500-50</b>	100 Pack

### Applications:

- General laboratory filtration
- Filtration for most HPLC samples

### Target2 PTFE Syringe Filters

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
4	0.2	No	<b>F2504-4</b>	100 Pack
	0.45	No	<b>F2504-3</b>	100 Pack
17	0.2	No	<b>F2513-4</b>	100 Pack
	0.45	No	<b>F2513-3</b>	100 Pack
30	0.2	No	<b>F2500-4</b>	100 Pack
		No	<b>F2500-3</b>	100 Pack
	Yes	<b>F2502-3</b>	100 Pack	
	1.0	No	<b>F2500-13</b>	100 Pack

### Applications:

- Filtration of aggressive organic, highly basic or hot solutions, transducer protectors
- Filter aqueous solutions after prewetting with an alcohol

### Target2 PVDF Syringe Filters

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
4	0.2	No	<b>F2504-6</b>	100 Pack
	0.45	No	<b>F2504-5</b>	100 Pack
17	0.2	No	<b>F2513-6</b>	100 Pack
	0.45	No	<b>F2513-5</b>	100 Pack
30	0.2	No	<b>F2500-6</b>	100 Pack
	0.45	No	<b>F2500-5</b>	100 Pack

### Applications:

- General biological filtration
- Filtration of samples where high protein recovery is desired



**Target2 Regenerated Cellulose Syringe Filters**

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
4	0.2	No	<b>F2504-8</b>	100 Pack
	0.45	No	<b>F2504-7</b>	100 Pack
17	0.2	No	<b>F2513-8</b>	100 Pack
	0.45	No	<b>F2513-7</b>	100 Pack
30	0.2	No	<b>F2500-8</b>	100 Pack
	0.45	No	<b>F2500-7</b>	100 Pack

**Target2 PES (Polyethersulfone) Syringe Filters**

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
17	0.2	No	<b>F2513-17</b>	100 Pack
	0.45	No	<b>F2513-14</b>	100 Pack
30	0.2	No	<b>F2500-17</b>	100 Pack
	0.45	No	<b>F2500-14</b>	100 Pack

**Target2 GMF (Glass Microfiber) Syringe Filters**

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
30	0.7	No	<b>F2500-18</b>	100 Pack
	1.2	No	<b>F2500-19</b>	100 Pack
	3.1	No	<b>F2500-20</b>	100 Pack

**Target2 Polypropylene Syringe Filters**

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
4	0.2	No	<b>F2504-10</b>	100 Pack
	0.45	No	<b>F2504-9</b>	100 Pack
17	0.2	No	<b>F2513-10</b>	100 Pack
	0.45	No	<b>F2513-9</b>	100 Pack
30	0.2	No	<b>F2500-10</b>	100 Pack
	0.45	No	<b>F2500-9</b>	100 Pack
		Yes	<b>F2502-9</b>	100 Pack

**Target2 Cellulose Acetate Syringe Filters**

Diameter (mm)	Pore Size (µm)	PreFilter	Cat. No.	Quantity
4	0.2	No	<b>F2504-16</b>	100 Pack
	0.45	No	<b>F2504-15</b>	100 Pack
17	0.2	No	<b>F2513-16</b>	100 Pack
	0.45	No	<b>F2513-15</b>	100 Pack
30	0.2	No	<b>F2500-16</b>	100 Pack
	0.45	No	<b>F2500-15</b>	100 Pack

**Applications:**

- Low non-specific binding applications
- Tissue culture media filtration and general biological sample filtration

**Applications:**

- Ion chromatography
- Tissue culture filtration, filtration of proteins and nucleic acids
- High-temperature liquids

**Applications:**

- Dissolution testing
- General filtration

**Applications:**

- Filtration of biological samples
- Filtration of aggressive organic solutions

**Applications:**

- Tissue culture media filtration, sensitive biological samples

## All-Plastic Disposable Syringes

Disposable syringes with polyethylene barrels and polypropylene plungers; use for all syringe filter applications

- Two-part, all-plastic construction eliminates the need for rubber or synthetic plunger gaskets
- No silicone or oil lubricant is required in the barrel
- Choose Luer-Slip or Luer Lock syringes, in capacities ranging from 1 to 50mL



### All-Plastic Disposable Syringes

Capacity (mL)	Cat. No.	Quantity
<b>Luer-Slip Syringes</b>		
1	<b>S7510-1</b>	100 Pack
3	<b>S7510-3</b>	100 Pack
5	<b>S7510-5</b>	100 Pack
10	<b>S7510-10</b>	100 Pack
20	<b>S7510-20</b>	100 Pack
30	<b>S7510-30</b>	50 Pack
50	<b>S7510-50</b>	30 Pack
<b>Luer Lock Syringes</b>		
3	<b>S7515-3</b>	100 Pack
5	<b>S7515-5</b>	100 Pack
10	<b>S7515-10</b>	100 Pack
20	<b>S7515-20</b>	100 Pack

## 750µL Micro-Centrifugal Filters, Non-Sterile

Filter volumes as low as 50µL with low hold-up volume

- Filter volumes as low as 50µL up to 750µL with low hold-up volume
- Use with any laboratory microcentrifuge
- Virgin polypropylene filter housing with tapered 2mL, capped receiver tube
- 10,000xG maximum centrifugal force



### 750µL Micro-Centrifugal Filters, Non-Sterile

Material	Pore Size (µm)	Cat. No.	Quantity
Cellulose Acetate	0.2	<b>F2517-1</b>	100 Pack
Cellulose Acetate	0.45	<b>F2517-2</b>	100 Pack
Nylon	0.2	<b>F2517-3</b>	100 Pack
Nylon	0.45	<b>F2517-4</b>	100 Pack
PVDF	0.2	<b>F2517-5</b>	100 Pack
PVDF	0.45	<b>F2517-6</b>	100 Pack
Regenerated Cellulose	0.2	<b>F2517-7</b>	100 Pack
Regenerated Cellulose	0.45	<b>F2517-8</b>	100 Pack
PTFE	0.2	<b>F2517-9</b>	100 Pack
PTFE	0.45	<b>F2517-10</b>	100 Pack

## 2mL Centrifugal Filters, Non-Sterile

- Filter sample volumes up to 2mL
- Virgin polypropylene filter housing with tapered 5mL, capped receiver tube
- Use with benchtop or floor model centrifuges
- 5,000xG maximum centrifugal force



### 2mL Centrifugal Filters, Non-Sterile

Material	Pore Size (µm)	Cat. No.	Quantity
Cellulose Acetate	0.2	<b>F2520-1</b>	25 Pack
Cellulose Acetate	0.45	<b>F2520-2</b>	25 Pack
Nylon	0.2	<b>F2520-3</b>	25 Pack
Nylon	0.45	<b>F2520-4</b>	25 Pack
PVDF	0.2	<b>F2520-5</b>	25 Pack
PVDF	0.45	<b>F2520-6</b>	25 Pack
PTFE	0.2	<b>F2520-7</b>	25 Pack
PTFE	0.45	<b>F2520-8</b>	25 Pack

## 25mL Centrifugal Filters, Non-Sterile

- Filter sample volumes up to 25mL
- Virgin polypropylene filter housing with conical receiver
- Use with benchtop or floor model centrifuges
- 2,500xG maximum centrifugal force



### 25mL Centrifugal Filters, Non-Sterile

Material	Pore Size (µm)	Cat. No.	Quantity
Cellulose Acetate	0.2	<b>F2519-1</b>	50 Pack
Cellulose Acetate	0.45	<b>F2519-2</b>	50 Pack
Nylon	0.2	<b>F2519-3</b>	50 Pack
Nylon	0.45	<b>F2519-4</b>	50 Pack
PVDF	0.2	<b>F2519-5</b>	50 Pack
PVDF	0.45	<b>F2519-6</b>	50 Pack

# Terms and Conditions of Sale

UNLESS OTHERWISE EXPRESSLY AGREED IN WRITING, ALL SALES ARE SUBJECT TO THE FOLLOWING TERMS AND CONDITIONS:

1. **GENERAL.** Thermo Fisher Scientific ("Seller") hereby offers for sale to the buyer named on the face hereof ("Buyer") the products listed on the face hereof (the "Products") on the express condition that Buyer agrees to accept and be bound by the terms and conditions set forth herein. Any provisions contained in any document issued by Buyer are expressly rejected and if the terms and conditions in this Agreement differ from the terms of Buyer's offer, this document shall be construed as a counter offer and shall not be effective as an acceptance of Buyer's document. Buyer's receipt of Products or Seller's commencement of the services provided hereunder will constitute Buyer's acceptance of this Agreement. This is the complete and exclusive statement of the contract between Seller and Buyer with respect to Buyer's purchase of the Products. No waiver, consent, modification, amendment or change of the terms contained herein shall be binding unless in writing and signed by Seller and Buyer. Seller's failure to object to terms contained in any subsequent communication from Buyer will not be a waiver or modification of the terms set forth herein. All orders are subject to acceptance in writing by an authorized representative of Seller.

2. **PRICE.** All prices published by Seller or quoted by Seller's representatives may be changed at any time without notice. All prices quoted by Seller or Seller's representatives are valid for thirty (30) days, unless otherwise stated in writing. All prices for the Products will be as specified by Seller or, if no price has been specified or quoted, will be Seller's price in effect at the time of shipment. All prices are subject to adjustment on account of specifications, quantities, raw materials, cost of production, shipment arrangements or other terms or conditions which are not part of Seller's original price quotation.

3. **TAXES AND OTHER CHARGES.** Prices for the Products exclude all sales, value added and other taxes and duties imposed with respect to the sale, delivery, or use of any Products covered hereby, all of which taxes and duties must be paid by Buyer. If Buyer claims any exemption, Buyer must provide a valid, signed certificate or letter of exemption for each respective jurisdiction.

4. **TERMS OF PAYMENT.** Seller may invoice Buyer upon shipment for the price and all other charges payable by Buyer in accordance with the terms on the face hereof. If no payment terms are stated on the face hereof, payment shall be net thirty (30) days from the date of invoice. If Buyer fails to pay any amounts when due, Buyer shall pay Seller interest thereon at a periodic rate of one and one-half percent (1.5%) per month (or, if lower, the highest rate permitted by law), together with all costs and expenses (including without limitation reasonable attorneys' fees and disbursements and court costs) incurred by Seller in collecting such overdue amounts or otherwise enforcing Seller's rights hereunder. Seller reserves the right to require from Buyer full or partial payment in advance, or other security that is satisfactory to Seller, at any time that Seller believes in good faith that Buyer's financial condition does not justify the terms of payment specified. All payments shall be made in U.S. Dollars.

5. **DELIVERY; CANCELLATION OR CHANGES BY BUYER.** The Products will be shipped to the destination specified by Buyer, F.O.B. Seller's shipping point. Seller will have the right, at its election, to make partial shipments of the Products and to invoice each shipment separately. Seller reserves the right to stop delivery of Products in transit and to withhold shipments in whole or in part if Buyer fails to make any payment to Seller when due or otherwise fails to perform its obligations hereunder. All shipping dates are approximate only, and Seller will not be liable for any loss or damage resulting from any delay in delivery or failure to deliver which is due to any cause beyond Seller's reasonable control. In the event of a delay due to any cause beyond Seller's reasonable control, Seller reserves the right to terminate the order or to reschedule the shipment within a reasonable period of time, and Buyer will not be entitled to refuse delivery or otherwise be relieved of any obligations as the result of such delay. Products as to which delivery is delayed due to any cause within Buyer's control may be placed in storage by Seller at Buyer's risk and expense and for Buyer's account. Orders in process may be canceled only with Seller's written consent and upon payment of Seller's cancellation charges. Orders in process may not be changed except with Seller's written consent and upon agreement by the parties as to an appropriate adjustment in the purchase price therefor. Credit will not be allowed for Products returned without the prior written consent of Seller.

6. **TITLE AND RISK OF LOSS.** Notwithstanding the trade terms indicated above and subject to Seller's right to stop delivery of Products in transit, title to and risk of loss of the Products will pass to Buyer upon delivery of possession of the Products by Seller to the carrier; provided, however, that title to any software incorporated within or forming a part of the Products shall at all times remain with Seller or the licensor(s) thereof, as the case may be.

7. **WARRANTY.** Seller warrants that the Products will operate or perform substantially in conformance with Seller's published specifications and be free from defects in material and workmanship, when subjected to normal, proper and intended usage by properly trained personnel, for the period of time set forth in the product documentation, published specifications or package inserts. If a period of time is not specified in Seller's product documentation, published specifications or package inserts, the warranty period shall be one (1) year from the date of shipment to Buyer for equipment and ninety (90) days for all other products (the "Warranty Period"). Seller agrees during the Warranty Period, to repair or replace, at Seller's option, defective Products so as to cause the same to operate in substantial conformance with said published specifications; provided that Buyer shall (a) promptly notify Seller in writing upon the discovery of any defect, which notice shall include the product model and serial number (if applicable) and details of the warranty claim; and (b) after Seller's review, Seller will provide Buyer with service data and/or a Return Material Authorization ("RMA"), which may include biohazard decontamination procedures and other product-specific handling instructions, then, if applicable, Buyer may return the defective Products to Seller with all costs prepaid by Buyer. Replacement parts may be new or refurbished, at the election of Seller. All replaced parts shall become the property of Seller. Shipment to Buyer of repaired or replacement Products shall be made in accordance with the Delivery provisions of the Seller's Terms and Conditions of Sale. Consumables are expressly excluded from this warranty.

Notwithstanding the foregoing, Products supplied by Seller that are obtained by Seller from an original manufacturer or third party supplier are not warranted by Seller, but Seller agrees to assign to Buyer any warranty rights in such Product that Seller may have from the original manufacturer or third party supplier, to the extent such assignment is allowed by such original manufacturer or third party supplier.

In no event shall Seller have any obligation to make repairs, replacements or corrections required, in whole or in part, as the result of (i) normal wear and tear, (ii) accident, disaster or event of force majeure, (iii) misuse, fault or negligence of or by Buyer, (iv) use of the Products in a manner for which they were not designed, (v) causes external to the Products such as, but not limited to, power failure or electrical power surges, (vi) improper storage and handling of the Products or (vii) use of the Products in combination with equipment or software not supplied by Seller. If Seller determines that Products for which Buyer has requested warranty services are not covered by the warranty hereunder, Buyer shall pay or reimburse Seller for all costs of investigating and responding to such request at Seller's then prevailing time and materials rates. If Seller provides repair services or replacement parts that are not covered by this warranty, Buyer shall pay Seller therefor at Seller's then prevailing time and materials rates. ANY INSTALLATION, MAINTENANCE, REPAIR, SERVICE, RELOCATION OR ALTERATION TO OR OF, OR OTHER TAMPERING WITH, THE PRODUCTS PERFORMED BY ANY PERSON OR ENTITY OTHER THAN SELLER WITHOUT SELLER'S PRIOR WRITTEN APPROVAL, OR ANY USE OF REPLACEMENT PARTS NOT SUPPLIED BY SELLER, SHALL IMMEDIATELY VOID AND CANCEL ALL WARRANTIES WITH RESPECT TO THE AFFECTED PRODUCTS.

THE OBLIGATIONS CREATED BY THIS WARRANTY STATEMENT TO REPAIR OR REPLACE A DEFECTIVE PRODUCT SHALL BE THE SOLE REMEDY OF BUYER IN THE EVENT OF A DEFECTIVE PRODUCT. EXCEPT AS EXPRESSLY PROVIDED IN THIS WARRANTY STATEMENT, SELLER DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, ORAL OR WRITTEN, WITH RESPECT TO THE PRODUCTS, INCLUDING WITHOUT LIMITATION ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. SELLER DOES NOT WARRANT THAT THE PRODUCTS ARE ERROR-FREE OR WILL ACCOMPLISH ANY PARTICULAR RESULT.

8. **INDEMNIFICATION.**

8.1 By Seller. Seller agrees to indemnify, defend and save Buyer, its officer, directors, and employees from and against any and all damages, liabilities, actions, causes of action, suits, claims, demands, losses, costs and expenses (including without limitation reasonable attorney's fees) ("Indemnified Items") for (i) injury to or death of persons or damage to property to the extent caused by the negligence or willful misconduct of Seller, its employees, agents or representatives or contractors in connection with the performance of services at Buyer's

premises under this Agreement and (ii) claims that a Product infringes any valid United States patent, copyright or trade secret; provided, however, Seller shall have no liability under this Section to the extent any such Indemnified Items are caused by either (i) the negligence or willful misconduct of Buyer, its employees, agents or representatives or contractors, (ii) by any third party, (iii) use of a Product in combination with equipment or software not supplied by Seller where the Product would not itself be infringing, (iv) compliance with Buyer's designs, specifications or instructions, (v) use of the Product in an application or environment for which it was not designed or (vi) modifications of the Product by anyone other than Seller without Seller's prior written approval. Buyer shall provide Seller prompt written notice of any third party claim covered by Seller's indemnification obligations hereunder. Seller shall have the right to assume exclusive control of the defense of such claim or, at the option of the Seller, to settle the same. Buyer agrees to cooperate reasonably with Seller in connection with the performance by Seller of its obligations in this Section.

Notwithstanding the above, Seller's infringement related indemnification obligations shall be extinguished and relieved if Seller, at its discretion and at its own expense (a) procures for Buyer the right, at no additional expense to Buyer, to continue using the Product; (b) replaces or modifies the Product so that it becomes non-infringing, provided the modification or replacement does not adversely affect the specifications of the Product; or (c) in the event (a) and (b) are not practical, refund to Buyer the amortized amounts paid by Buyer with respect thereto, based on a five (5) year amortization schedule. THE FOREGOING INDEMNIFICATION PROVISION STATES SELLER'S ENTIRE LIABILITY TO BUYER FOR THE CLAIMS DESCRIBED HEREIN.

8.2 By Buyer. Buyer shall indemnify, defend with competent and experienced counsel and hold harmless Seller, its parent, subsidiaries, affiliates and divisions, and their respective officers, directors, shareholders and employees, from and against any and all damages, liabilities, actions, causes of action, suits, claims, demands, losses, costs and expenses (including without limitation reasonable attorneys' fees and disbursements and court costs) to the extent arising from or in connection with (i) the negligence or willful misconduct of Buyer, its agents, employees, representatives or contractors; (ii) use of a Product in combination with equipment or software not supplied by Seller where the Product itself would not be infringing; (iii) Seller's compliance with designs, specifications or instructions supplied to Seller by Buyer; (iv) use of a Product in an application or environment for which it was not designed; or (v) modifications of a Product by anyone other than Seller without Seller's prior written approval.

9. **SOFTWARE.** With respect to any software products incorporated in or forming a part of the Products hereunder, Seller and Buyer intend and agree that such software products are being licensed and not sold, and that the words "purchase", "sell" or similar or derivative words are understood and agreed to mean "license", and that the word "buying" or similar or derivative words are understood and agreed to mean "licensing". Notwithstanding anything to the contrary contained herein, Seller or its licensor, as the case may be, retains all rights and interest in software products provided hereunder.

Seller hereby grants to Buyer a royalty-free, non-exclusive, nontransferable license, without power to sublicense, to use software provided hereunder solely for Buyer's own internal business purposes on the hardware products provided hereunder and to use the related documentation solely for Buyer's own internal business purposes. This license terminates when Buyer's lawful possession of the hardware products provided hereunder ceases, unless earlier terminated as provided herein. Buyer agrees to hold in confidence and not to sell, transfer, license, loan or otherwise make available in any form to third parties the software products and related documentation provided hereunder. Buyer may not disassemble, decompile or reverse engineer, copy, modify, enhance or otherwise change or supplement the software products provided hereunder without Seller's prior written consent. Seller will be entitled to terminate this license if Buyer fails to comply with any term or condition herein. Buyer agrees, upon termination of this license, immediately to return to Seller all software products and related documentation provided hereunder and all copies and portions thereof.

Certain of the software products provided by Seller may be owned by one or more third parties and licensed to Seller. Accordingly, Seller and Buyer agree that such third parties retain ownership of and title to such software products. The warranty and indemnification provisions set forth herein shall not apply to software products owned by third parties and provided hereunder.

10. **LIMITATION OF LIABILITY.** NOTWITHSTANDING ANYTHING TO THE CONTRARY CONTAINED HEREIN, THE LIABILITY OF SELLER UNDER THESE TERMS AND CONDITIONS (WHETHER BY REASON OF BREACH OF CONTRACT, TORT, INDEMNIFICATION, OR OTHERWISE, BUT EXCLUDING LIABILITY OF SELLER FOR BREACH OF WARRANTY (THE SOLE REMEDY FOR WHICH SHALL BE AS PROVIDED UNDER SECTION 7 ABOVE)) SHALL NOT EXCEED AN AMOUNT EQUAL TO THE LESSER OF (A) THE TOTAL PURCHASE PRICE THEREFORE PAID BY BUYER TO SELLER WITH RESPECT TO THE PRODUCT(S) GIVING RISE TO SUCH LIABILITY OR (B) ONE MILLION DOLLARS (\$1,000,000). NOTWITHSTANDING ANYTHING TO THE CONTRARY CONTAINED HEREIN, IN NO EVENT SHALL SELLER BE LIABLE FOR ANY INDIRECT, SPECIAL, CONSEQUENTIAL OR INCIDENTAL DAMAGES (INCLUDING WITHOUT LIMITATION DAMAGES FOR LOSS OF USE OF FACILITIES OR EQUIPMENT, LOSS OF REVENUE, LOSS OF DATA, LOSS OF PROFITS OR LOSS OF GOODWILL), REGARDLESS OF WHETHER SELLER (a) HAS BEEN INFORMED OF THE POSSIBILITY OF SUCH DAMAGES OR (b) IS NEGLIGENT.

11. **EXPORT RESTRICTIONS.** Buyer acknowledges that each Product and any related software and technology, including technical information supplied by Seller or contained in documents (collectively "Items"), is subject to export controls of the U.S. government. The export controls may include, but are not limited to, those of the Export Administration Regulations of the U.S. Department of Commerce (the "EAR"), which may restrict or require licenses for the export of Items from the United States and their re-export from other countries. Buyer shall comply with the EAR and all other applicable laws, regulations, laws, treaties, and agreements relating to the export, re-export, and import of any Item. Buyer shall not, without first obtaining the required license to do so from the appropriate U.S. government agency, (i) export or re-export any Item, or (ii) export, re-export, distribute or supply any Item to any restricted or embargoed country or to a person or entity whose privilege to participate in exports has been denied or restricted by the U.S. government. Buyer shall, if requested by Seller, provide information on the end user and end use of any Item exported by the Buyer or to be exported by the Buyer. Buyer shall cooperate fully with Seller in any official or unofficial audit or inspection related to applicable export or import control laws or regulations, and shall indemnify and hold Seller harmless from, in connection with, any violation of this Section by Buyer or its employees, consultants, or agents.

12. **MISCELLANEOUS.** (a) Buyer may not delegate any duties nor assign any rights or claims hereunder without Seller's prior written consent, and any such attempted delegation or assignment shall be void. (b) The rights and obligations of the parties hereunder shall be governed by and construed in accordance with the laws of the State of Seller's manufacturing location, without reference to its choice of law provisions. Each party hereby irrevocably consents to the exclusive jurisdiction of the state and federal courts located in the county and state of Seller's manufacturing location, in any action arising out of or relating to this Agreement and waives any other venue to which it may be entitled by domicile or otherwise. (c) In the event of any legal proceeding between the Seller and Buyer relating to this Agreement, neither party may claim the right to a trial by jury, and both parties waive any right they may have under applicable law or otherwise to a right to a trial by jury. Any action arising under this Agreement must be brought within one (1) year from the date that the cause of action arose. (d) The application to this Agreement of the U.N. Convention on Contracts for the International Sale of Goods is hereby expressly excluded. (e) In the event that any one or more provisions contained herein shall be held by a court of competent jurisdiction to be invalid, illegal or unenforceable in any respect, the validity, legality and enforceability of the remaining provisions contained herein shall remain in full force and effect, unless the reformer materially changes the bargain. (f) Seller's failure to enforce, or Seller's waiver of a breach of, any provision contained herein shall not constitute a waiver of any other breach or of such provision. (g) Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals. (h) Buyer agrees that all pricing, discounts and technical information that Seller provides to Buyer are the confidential and proprietary information of Seller. Buyer agrees to (1) keep such information confidential and not disclose such information to any third party, and (2) use such information solely for Buyer's internal purposes and in connection with the Products supplied hereunder. Nothing herein shall restrict the use of information available to the general public. (i) Any notice or communication required or permitted hereunder shall be in writing and shall be deemed received when personally delivered or three (3) business days after being sent by certified mail, postage prepaid, to a party at the address specified herein or at such other address as either party may from time to time designate to the other.

# Resources

## for Chromatographers



### Chromatography Resource Center

Our web-based resource center provides technical support, applications, technical tips and literature to help move your separations forward.

Visit [www.thermoscientific.com/chromatography](http://www.thermoscientific.com/chromatography)

### How to order

#### USA and Canada

800 332 3331  
865 354 4616 fax  
[chromatography@thermofisher.com](mailto:chromatography@thermofisher.com)

#### France

+33 (0)1 60 92 48 34  
+33 (0)1 60 92 49 00 fax  
[Consommableslscgcf@thermofisher.com](mailto:Consommableslscgcf@thermofisher.com)

#### Germany

+49 (0) 2423 9431 -20 or -21  
+49 (0) 6103 408 1640 fax  
[lcgc.consumables.DE@thermofisher.com](mailto:lcgc.consumables.DE@thermofisher.com)

#### United Kingdom

+44 1928 534110  
+44 1928 534001 fax  
[lcgc.consumables@thermofisher.com](mailto:lcgc.consumables@thermofisher.com)

#### China

+86 21 6865 4588

#### India

+91 22 6742 9494  
+91 27 1766 2352  
+91 22 6742 9495 fax  
[analyze.in@thermofisher.com](mailto:analyze.in@thermofisher.com)

#### Japan

+81 3 5826 1615  
+81 3 5826 1611 fax  
[analyze.jp@thermofisher.com](mailto:analyze.jp@thermofisher.com)

#### Thermo Fisher Scientific Australia Pty Ltd

1300 735 292 (free call domestic)  
1800 067 639 fax  
[AUinfo@thermofisher.com](mailto:AUinfo@thermofisher.com)

#### Thermo Fisher Scientific New Zealand Ltd

0800 933 966 (free call domestic)  
0800 329 246 fax  
[NZinfo@thermofisher.com](mailto:NZinfo@thermofisher.com)

#### All Other Enquiries

+44 (0) 1928 534 050  
+44 (0) 1928 534 049 fax  
[salesorders.columns.uk@thermofisher.com](mailto:salesorders.columns.uk@thermofisher.com)

### Technical Support

**For advice and support,  
please visit our website:**

[www.thermoscientific.com/chromexpert](http://www.thermoscientific.com/chromexpert)

For more information visit: [www.thermoscientific.com/chromatography](http://www.thermoscientific.com/chromatography)

© 2015 Thermo Fisher Scientific Inc. All rights reserved.  
Specifications, terms and pricing are subject to change. Not all products are available in all countries.  
Please consult your local sales representative for details.

**Thermo**  
SCIENTIFIC

A Thermo Fisher Scientific Brand